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E-posta: info@htcert.com

Symbols and Their Explanations

-  Read the user manual before use!
-  Warning!
-  It is not sterile!
-  Do not reuse !
-  Store in a dry place!
-  Protect from sunlight!
-  Humidity limitation symbol: Store in environments not exceeding 0-80% humidity!
-  Manufacturer
-  CE marking according to (EU) 2017/745 Medical Device Directive

References

(EU) 2017/745 Medical Device Regulation

EN ISO 20417:2021 Medical devices - Information to be provided by the manufacturer

EN ISO 15223-1 Medical devices - Symbols to be used in the information to be submitted by the manufacturer - Part 1: General requirements

EN ISO 17664-1 Processing of health care products - Information provided by medical device manufacturers for the processing of medical devices - Part 1: Critical and semi-critical medical devices

EN ISO 17665-1 Sterilization of healthcare products - Moist heat - Part 1 - Requirements for the development, validation and routine control of a sterilization process for medical devices

EN ISO 17665-2 Sterilization of health care products - Moist heat – Part 2 - Requirements for the development, validation and routine control of a sterilization process for medical devices

EN ISO 14630:2012 Non-active surgical implants - General specifications

ISO 11607-1:2019 (Packaging for Terminally Sterilized Medical Devices - Part 1: Requirements for Materials, Sterile Barrier Systems and Packaging Systems)

Product Description

Ortimplant Cannulated Bone Screw Implants are manufactured as threaded for use alone or with other implant products for fracture fixation or correction of bone disorder. Screw implants are made of stainless steel conforming to ISO 5832-1 or titanium alloy conforming to ISO 5832-2 and are available in various sizes and designs. Titanium and stainless steel are completely biocompatible, corrosion resistant and non-toxic.

Ortimplant implantable screws are sold non-sterile and all devices are disposable.

Cannulated bone screws and components should only be used for bone fixation by qualified personnel with the necessary application competencies and full experience and certification in the medical field. These products, components, treatment kits and tools are manufactured for the purpose of medical interventions for the treatment of patients and may not be used for any other purposes. Otherwise, Ortimplant cannot be held responsible.

Note: The size, type and dimensions of the implant and accessories to be used, "application sets" suitable for the selected implant and accessories should be determined by the surgeon who will perform the application according to the bone area on which he will use, the condition of the disease, the gender of the patient and the age of the patient.

Before using Ortimplant products, surgeons and their assistant team should read these statements carefully. All necessary documents and materials must be provided by the authorized distributor. Accordingly, the surgeon should know the risks that may arise as a result of the use of these products and inform the patient about these risks as much as possible. Product implantation should be performed by experts trained and experienced in medical plate application.

Cannulated screws, cannulated screw system are generally designed for fracture fixation of bones. While it is used for plate fixation in fracture fixations, it can also be used alone as a free fragment or for reduction and compression in the cancellous area. It can be used for femoral neck fractures, intercondylar femoral, malleolus, pilon, tibia, calcaneus, talus, tibial plateau, tarsal, metatarsal, wrist, metacarpal, carpal, scaphoid and radius fractures and arthrodesis. It can also be used in the correction of iliosacral dislocations and hallux valgus.

Screw Washers are used to increase the bone contact area of the screw in cannulated screws used alone and to prevent the screw from being embedded in the bone when used at an angle.

Headless screws are used for adult and pediatric patients in the treatment of small fragment fractures (hand and foot bones, carpal, metacarpal, tarsal, metatarsal, phalangeal bones), distal and proximal humerus, radius and ulna fractures, osteotomy, nonunion, malunion, pseudoarthroses, arthrodeses. is used. It is used for temporary fixation of large fragment fractures (femoral neck fractures, intercondylar femur fractures, tibial plateau fractures), femoral head epiphysiolysis, ankle arthrodesis, and iliosacral dislocations. It is used in ligamentous injuries of small, medium and large joints.

Herbert screws are used for adult and pediatric patients in the treatment of hand and foot bones (carpal, metacarpal, tarsal, metatarsal, phalangeal bones), distal and proximal humerus, radius and ulna fractures, osteotomy, nonunion, malunion, pseudoarthroses, arthrodeses. Small medium It is also used in ligamentous injuries of large joints. It is especially used in the treatment of scaphoid bone fractures and pseudoarthroses.

Interference screws are used for temporary fixation of bone patellar tendon bone and soft tissue grafts in anterior cruciate ligament and posterior cruciate ligament reconstructions.

Malleolar screws are used in fractures, fixation of metatarsal fractures, malunions and unions, acute fractures, malleolar fractures, talus fractures, foot and ankle osteotomy in order to provide osteotomy and bone reconstruction.

Bone screw implants produced by Ortimplant have different diameters, lengths, profiles and tooth feed, and not every screw can be used with all plates, intramedullary nails, external fixator implants and components. Suitable screws and components (application sets) for the appropriate implant should be determined by the surgeon during preoperative planning. The materials should be sent for sterilization after the necessary compatibility checks are made by the surgeon and/or his team before the surgery. Otherwise, the surgeon is responsible for the surgery and postoperative complications. Ortimplant and its Distributors are not responsible to patients and surgeons under these circumstances.

Types of Cannulated Screw Implants and the surgical sets they are used with;

- ✓ 3.5 mm / 4.0 mm Cannulated Screws; Must be used with 3.5 / 4.0 mm Cannulated Screw Application Tool Sets,
- ✓ 6.5 / 7.3 mm Cannulated Screws; It should be used with 6.5 / 7.3 mm Cannulated Screw Application Tool Sets.
- ✓ Herbert Cannulated Screws; It should be used with Herbert Cannulated Screw Application Tool Sets.
- ✓ 3.5 / 6.5 mm Headless Cannulated Screw should be used with 3.5 / 6.5 mm Headless Cannulated Screw Application Tool Sets.

Specialists are responsible for the consequences of using the wrong implant, not meeting the safety criteria in this user guide, or using inappropriate surgical techniques. In such cases, Ortimplant and its authorized distributor are not responsible to the patient and surgeon.

Some types of Ortimplant cannulated bone screw implants can be used with implants and application sets produced by other companies, but this does not mean that all screws can be used. The surgeon should definitely avoid the use of such inappropriate products as he will not be able to provide solid bone fixation as a result of such problems that may arise.

During the preparation phase before the surgery, the surgeon must ensure that the screw implants, application sets and implants he will use are compatible with each other. Ortimplant and its distributors are not responsible for the negative consequences of surgical application in case of any complications.

 **Note:** Cannulated Bone Screws are manufactured from ISO 5832-1 stainless steel and ISO 5832-3 titanium alloy certified raw materials. Do not use titanium alloy and stainless steel screws in the same place.

 **Note:** To achieve the best result; Do not use Ortimplant cannulated screw implants with parts produced by other companies. Do not reuse bone screw implants due to metal fatigue and risk of infection.

 **Note to Surgeon:** Even if the physician is a knowledgeable intermediary between the patient and the company, the important medical information provided in this document must be disclosed to the patient.

Intended Use

All cannulated screw implants; They are used in the treatment of bone fractures in the body, congenital or acquired bone-joint disorders, depending on the type, model, size and surgeon's choice.

Indication Information

Cannulated Screws, Headless screws, Herbert screws, Interference screws, Malleolar screws and Washers (Washer) Ortimplant Cannulated Bone Screw Implants are designed to provide fixation for fractures, fusions or osteotomies. In general, these screws are indicated for fractures that require additional stability.

Indications include:

Cannulated Screws, Headless screws, Herbert screws, Interference screws, Malleolar screws and Washers (Washer); fracture of the hands, feet, pelvis and sacrum with metaphyseal fracture of the extremities

- ✓ 3.5 / 4.0 mm Cannulated Screws; It is used for fixation of small bone fractures.
- ✓ 6.5 / 7.3 mm Cannulated Screws; It is used for fixation of large bone fractures.
- ✓ Herbert Cannulated Screws; It is used for fixation of scaphoid fractures and small bone fractures.

Note: The size, type and dimensions of the implants and accessories to be used in surgery and the "application sets" suitable for the selected implants and accessories should be determined by the surgeon who will perform the application according to the bone area to be used, the structure of the bone, the condition of the disease, the gender of the patient and the age of the patient. The Ortimeplant catalog contains information about the plates to be used with implants.

Planning the surgical procedure should be based on detailed evaluation of the patient's x-rays. In this way, the appropriate bone screw implant selection is determined by the surgeon by calculating size and finish combinations. Pre-operative planning templates are available from the manufacturer. The patient's allergic sensitivity to the implant material should be determined. Implant screws of different sizes should be available at the time of surgery. Before placing the appropriate implant, the surgeon must ensure that the implant size and type are compared to the label to use the correct size implant. Excess labels inside the box; for reporting and ordering purposes.

Note: Cannulated Bone screw implants, components and application sets suitable for the implant should be selected by the surgeon and his team before the operation, and their operation should be checked and verified.

The surgeon and his team are responsible for selecting the appropriate implant and its components for the patient and the application sets suitable for these implants and components. Ortimeplant and distributors cannot be held responsible for any problems arising from the selection and application of inappropriate implants, components and application sets.

Contraindication Information

The following situations include, but are not limited to, contraindications;

- ✓ Development of active infection or risk of infection (immunodeficiency),
- ✓ Indicators of local inflammation,
- ✓ Increased fever and white blood cell mass,
- ✓ Circulatory disorders and occlusions that prevent adequate blood circulation in the fracture or surgical area,
- ✓ Morbid obesity,
- ✓ Pregnancy (until the patient's suitability for surgery conditions is confirmed by experts),
- ✓ mental illnesses,
- ✓ Some mental or neuromuscular diseases that may cause complications or fixation failure in the postoperative period,
- ✓ Patients with serious congenital anatomical disorders,
- ✓ The presence of another disease or associated use of a medication or device may potentially reduce or prevent the function of screw implants (e.g., progressive streptococcal infections, congenital functional abnormality, unexplained semimetal elevation, other disabling conditions).
- ✓ Bone stock disorders that prevent fixation or support of the device,
- ✓ Patients with any anomalies in the bone anatomy in the area where the screw is intended to be implanted,
- ✓ Rapidly progressive joint disease, bone erosion, osteopenia and/or osteoporosis. Osteoporosis and osteopenia are contraindications that may limit the degree of healing due to poor stabilization and/or mechanical fixation.
- ✓ Cases of metal allergy to titanium and stainless steel,
- ✓ Situations where bone graft or fusion is unnecessary,
- ✓ In cases where the implant components selected for use are too large or too small for the bone area to be applied,
- ✓ Patients who do not have sufficient soft tissue, sufficient bone stock, or sufficient bone quality to cover the operation area,
- ✓ Patients expected to have low physiological performance,
- ✓ Patients who are not willing to follow postoperative instructions
- ✓ Using the products in situations not specified in the indications,

Note: Due to the high risk of unexpected metal corrosion inside the body, two different metal components should not be used together.

Potential adverse events and complications

All available scientific safety warnings regarding bone surgery should be taken into account by the surgeon and the team during the placement of the bone screw implant. In many cases, negative results may be clinic-related rather than device-related.

Potential safety warnings, potential adverse events, and complications related to implantation of bone screw implants and components are listed below, but are not limited to;

- ✓ Nonunion or delayed union of the fracture,
- ✓ Early or late loosening of bone screw implants and components,
- ✓ Possibility of disassembly, bending and/or breakage of all components within the body,
- ✓ Development of foreign body reaction (allergy) to implants, rash formation, corrosion of implants (cracks, friction and/or general corrosion)
- ✓ In patients without adequate soft tissue cover, skin penetration due to pressure on the operation side, increase in fibrous tissue, development of skin necrosis and/or pain, and inflammation of surrounding tissues and bone may occur. Improper use of instruments or incorrect placement of the implant may cause damage to surrounding tissue and nerves.
- ✓ After the surgery, changes in bone smoothness, limb weakening and/or length loss in the relevant body part should be monitored by the surgeon and the patient.
- ✓ Development of local and/or general infection,
- ✓ Paralysis (complete or partial), dysesthesia, hypersensitivity, numbness, paresthesia, appearance of radiculopathy and/or pain, depersonalization, loss of neurological function such as neuroma development, loss of consciousness, tingling sensation and/or visible defects,
- ✓ Disorders of bladder fullness, urinary retention and other urinary system disorders,
- ✓ Scar formation can cause nervous system problems or pressure and pain around the nerve,
- ✓ Fracture, microfracture, resorption, development of bone defect or penetration of bone,
- ✓ Pseudoarthrosis, late and poor union of the fractured bone where the surgical operation was performed,
- ✓ Increased/decreased mobility and function in the adjacent joints of the operated bone,
- ✓ The patient cannot continue his daily activities,
- ✓ Loss or decrease in bone density (possibly a result of a decrease in the load on the bone),
- ✓ Complications of the grafted area, including pain, fracture or wound healing problems,
- ✓ Intestinal obstruction or loss of control, gastritis or other gastrointestinal tract disorders
- ✓ deep vein thrombosis,
- ✓ Bleeding, hematoma, seroma, edema, hypertension, vascular occlusion, stroke, excessive bleeding, phlebitis, tissue decay, wound dehiscence, vascular injuries or other cardiovascular system disorders,
- ✓ Avascular necrosis,
- ✓ Shortness in the limb where the implant was applied,
- ✓ Reproductive system disorders, sexual dysfunction, progression of respiratory problems (e.g. pulmonary embolism, pneumonia, bronchitis),
- ✓ Damage to the epiphysis of long bones as a result of placing inappropriate products, wrong size screws and guide wires during surgery,
- ✓ Patients should stay away from strong magnetic fields.
- ✓ The patient should immediately consult his physician if he hears abnormal sounds, pain, or increased temperature in the surgery area.
- ✓ It has been reported that during implantation of screws, especially around joints, inappropriately sized screws can penetrate joints and cause chondral lesions. It is necessary to perform the surgical procedure precisely to prevent chondral lesions and early joint arthrosis.

Note: Additional surgical intervention may be required to correct these potential disadvantages.

Intended User Profile

Ortimplant Cannulated Screws should only be used by healthcare professionals.

Surgical procedures should only be performed by medical specialists or medical residents under the responsibility of specialists who have received adequate training in the anatomical structure or formation where the surgery will be performed and are familiar with surgical techniques.

Before performing a surgical procedure, surgical literature on techniques, complications, and hazards should be reviewed. You should carefully read all operating instructions regarding product safety features before use. All necessary materials and documents must be provided by the authorized distributor.

General Safety Warnings

- ! Read the User Manual before use and take into account the warnings.
- ! The products are NOT STERILE! Before first use, sterilize according to the instructions specified in this manual.
- ! Do not use it for purposes other than its intended use.
- ! Implant cannulated screws should not be used for purposes other than those stated above.
- ! Choosing the appropriate device, appropriate accessory and appropriate application tool set for the patient is very important. Improper use or misalignment of components may cause loss, bending, breakage, or cracking of the device or bones.
- ! Until radiographic confirmation of adequate bone development; Other internal or external support devices should be used (orthosis, plaster, etc.) to stabilize the fracture and minimize the load on the fractured area.
- ! Fractures and osteotomies of long bones close to the joints increase the stress on the bone screws. A sufficient number of screws should be used in fracture fixation to ensure maximum fixation. The duration of weight bearing on the fracture site should be extended until adequate bone union occurs.
- ! Screw implants should not be placed in the fracture line. The screw threads should fit snugly into the bone and the screws should be long enough to allow telescopic movement. However, it should be adjusted so that it does not protrude too much from the bone cortex and cause vascular nerve injury or damage to the joint surfaces.
- ! Do not mix incompatible metals.
- ! The device is made of metal materials and does not contain waste materials.

Precautions

Implant cannulated screws should not be reused. Loading too early may cause devices to stop working. In order to prevent the implant from causing discomfort during use, check the implants and application sets (in terms of suitability and full functionality) before use. Precautions to be taken before, during and after other procedures are listed below.

Operation

Patient and Implant Selection:

Only patients with indications and no absolute contraindications should undergo surgical procedure. However, the final decision and responsibility belongs to the surgeon who will perform the surgery.

Choosing the appropriate size, shape and design for each patient is important for the procedure to be successful. Metallic surgical implants are subjected to repeated loads during use and are limited in shape, size and design. Unless care is taken in selecting the appropriate implant for the patient, unless the implant is placed properly, and if no precautions are taken to reduce the stress on the implant, such stresses can result in buckling and device loosening, metal fatigue, and eventual fracture before healing. These situations may result in device damage or premature removal.

Before Operation:

- ✓ Only patients who meet the criteria specified in the indications should be selected.
- ✓ The use of screw implants should be avoided in patients who have or are prone to the conditions specified in the contraindications.
- ✓ Care should be taken to transport and store implant components. Screw implants should not be scratched or damaged. Screws must be protected during storage against aggressive environmental conditions. Screw implant packages must be undamaged.
- ✓ Screws of appropriate size and number should be used. Additionally, excess stock and variety should be kept ready during surgery.

- ✓ The surgeon must become familiar with the mechanical parts of the delivery sets, implants, and various other components before using the surgical equipment and must set up and check the system to verify that all parts of the implants, delivery sets, and other components, and necessary instrumentation are present and working harmoniously before the procedure begins. Bone screw implants cannot be completed with medical devices manufactured by other manufacturers. Different types of metal should never be used together.
- ✓ Components and instruments must be clean and sterile before use. Additional sterile materials (implant size, type alternatives, application sets, hospital orthopedic surgery sets, etc.) should be available to prevent undesirable situations.

During the operation:

- ✓ Extreme care should be taken in the surgical approach to the areas through which the spinal cord, nerve roots and peripheral nervous system pass. Damage to nerves leads to loss of neurological functions. During the surgical procedure, maximum attention should be paid to intracranial, intrathoracic, intrapelvic organs, arterial and venous areas and these structures should be protected. Arterial or venous injuries can result in death.
- ✓ Breakage, slippage, or misuse of instrument sets or implant components may result in injury to the patient or staff.
- ✓ Bone screw implants should not be bent repeatedly. Since scratches or nicks on the implant may reduce the functional integrity of the implant, due care should be taken by the surgeon to avoid use of the implant in such cases.
- ✓ Whenever possible or necessary, an imaging system should be used to facilitate surgery.
- ✓ Do not use bone bonding agents and solvents as their safety and effectiveness have not been proven as it will be difficult to separate the implant components from each other. Drills and cutters used during drilling and cutting operations cause high heat due to friction. Heat accumulation occurring during surgery can cause neurological injuries and local bone tissue death.
- ✓ Before suturing the soft tissue, change the screw implant position if there is an inappropriate screw implant position.

After the operation:

- ✓ The patient should be informed to avoid placing excessive regional load on the operation area. In such cases, before bone union is completed, the patient should be informed that bone screw implants may bend and bend, and the fracture position may be distorted. The patient should be informed by the surgeon about the risk of reoperation under these conditions.
- ✓ The surgeon should inform the patient about contraindications that may cause the device to loosen and/or break. If the patient is extremely elderly, frail, or mobile, the risk of bending, loosening, and breaking of the temporary fixation device may be increased during the postoperative rehabilitation period. The patient should be warned to avoid falls and sudden jolts in different positions.
- ✓ To achieve maximum surgical success, neither the patient nor the device should be subjected to mechanical vibration or shock that could cause the device to loosen. The patient should be informed about this possibility and, in particular, his/her participation in any sports activities should be restricted. The patient should be advised to be accompanied by those nearby and to use auxiliary support devices (crutches, walkers, orthoses, etc.) while standing up. The patient should especially be advised to avoid turning movements. The patient should be warned not to drink or consume a product containing nicotine, alcohol, or steroids during the bone healing process.
- ✓ The patient should be warned about the loss of movement in the surgical area and informed about how to compensate for the permanent physical limitation in body movement.
- ✓ Nonunion or delayed union of the fracture; creates excessive and repetitive stress on the implant that stabilizes the fracture. Excessive load may cause the metal device to bend, loosen, or collapse due to fatigue. It is important not to move the surgical area until appropriate bone union is achieved by x-ray examination. If the device is loose, bent and/or warped, it must be replaced or removed without serious injury. The patient should be warned of these dangers and monitored closely until the bone consolidates.
- ✓ Especially in high-risk patients, prophylactic antibiotic use is precautionary before subsequent surgical procedures (such as dental surgery).
- ✓ Implant fixation devices are designed to provide stability at the fracture site during normal fracture healing. Once the bone is fused, these devices have no functional effect and can be removed. The final decision to remove the device is up to the surgeon and the patient. Many patients may experience one or more of the following complications;

- Corrosion, local tissue reaction or pain,
- Disruption of implant position (leading to possible injuries),
- Extra risk of injury due to trauma to the surgical area after surgery,
- Bending, loosening or breaking of the implant in the surgical area, making it difficult to remove it,
- Pain, discomfort or abnormal movement while the implant is in the body,
- High risk of possible infection,
- Bone loss caused by excessive pressure on the bone.

Removing the Implant:

The length of time the device remains on the patient is determined by the surgeon, taking into account the progress of the union and the duration of treatment. Before the screws are removed from the implants, the patient is given an appropriate dose of paracetamol. Turning the bone screws counterclockwise easily loosens them and pulls them out of the bone. The wounds are closed by placing dry sterile gauze on them.

For detailed information about removing screws, please contact your doctor.

Note: Ortimplant; recommends removing the implant after the fracture has healed to prevent refracture of the bone and other complications. The bone screw implant and its components should not be reused. It is not possible to use a device removed from one patient on another patient and its use should not even be considered.

Surgical Technique

A separate surgical technique guide is provided to explain the procedure for implant screws during the surgical procedure. The surgeon should be familiar with the surgical procedure before using the products. It is also the surgeon's responsibility to consult with colleagues experienced in the procedures and to become familiar with relevant publications before use. Please review the surgical technique document provided by Ortimplant for the correct application method.

Sterilization

Steam Sterilization Method has been validated and is recommended as the sterilization method for our products.

- ! The steam sterilizer must comply with the EN 285 standard and be validated according to the EN ISO 17665 standard.
- ! It must be ensured that the maximum allowable loading is not exceeded according to the sterilizer manufacturer's instructions.
- ! The properties of the water to be used in the sterilization process must comply with EN 285 and ISO 17665 standards.

Sterilization Method	Temperature	Time
Steam Sterilization	134 °C	5 min

Packaging and Storage Conditions

All Cannulated Screws are cleaned and disinfected according to standards-compliant methods. Cleaning and decontamination involves rinsing with neutral cleaning agents followed by washing with distilled water.

Afterwards, all products are packaged in the clean room.

Note: It should be noted that the packages containing Ortimplant Cannulated Screws are not damaged. If a delivery note or invoice is used, check that all products are complete, including application sets, before use. Damaged packages and products should not be used and should be returned to Ortimplant or the distributor.

Cannulated Screws should be stored in an environment with humidity below 80% and the environment should not contain corrosive gas.

Disposal Conditions

Disposal of the products must be carried out in accordance with the national legislation of the country where the products are used.

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Manufacturer: Ortimplant Ortopedi ve Tıbbi Aletler Tic. Ltd. Şti.