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ENGLISH

Symbols and Their Explanations

-  Read the user manual before use!
-  Warning!
-  It is not sterile!
-  Do not reuse !
-  Store in a dry place!
-  Protect from sunlight!
-  Humidity limitation symbol: Store in environments not exceeding 0-80% humidity!
-  Manufacturer
-  **CE**₂₈₀₃ CE marking according to (EU) 2017/745 Medical Device Directive

References

(EU) 2017/745 Medical Device Regulation**EN ISO 20417:2021** Medical devices - Information to be provided by the manufacturer**EN ISO 15223-1** Medical devices - Symbols to be used in the information to be submitted by the manufacturer - Part 1: General requirements**EN ISO 17664-1** Processing of health care products - Information provided by medical device manufacturers for the processing of medical devices - Part 1: Critical and semi-critical medical devices**EN ISO 17665-1** Sterilization of healthcare products - Moist heat - Part 1 - Requirements for the development, validation and routine control of a sterilization process for medical devices**EN ISO 17665-2** Sterilization of health care products - Moist heat – Part 2 - Requirements for the development, validation and routine control of a sterilization process for medical devices**EN ISO 14630:2012** Non-active surgical implants - General specifications**ISO 11607-1:2019** (Packaging for Terminally Sterilized Medical Devices - Part 1: Requirements for Materials, Sterile Barrier Systems and Packaging Systems)**Product Description**

Ortimplant External Fixators should only be used by specialist healthcare professionals (surgeons) who have the necessary application qualifications for the treatment of orthopedic disorders and traumatology and have full experience and certificates in this field. These products, components, treatment kits and tools are manufactured for the purpose of medical interventions for the treatment of patients and may not be used for any other purposes. Otherwise, Ortimplant cannot be held responsible.

Ortimplant External Fixators and their components are sold non-sterile and all devices are disposable.

External Fixators are manufactured from certified raw materials such as aluminum, stainless steel, fiber carbon, titanium and titanium alloy. Ortimplant devices are manufactured from materials that provide the above features for the purpose. Ortimplant Bone External Fixators are manufactured from aluminum, titanium, fiber carbon and stainless steel materials. They are completely biocompatible, corrosion resistant and non-toxic in the biological environment.

Note: The size, type and dimensions of the external fixator implants and accessories to be used and the "application sets" suitable for the selected implant and accessories should be determined by the surgeon who will perform the application according to the bone area on which he will use, the condition of the disease, the gender of the patient and the age of the patient.

Before using the Ortimplant external fixator and its components, surgeons and the supporting team should carefully read the instructions for use. All necessary documents and materials will be provided to you by the authorized distributor. The surgeon should be aware of the risks that may arise from the use of these products and inform the patient about these risks as much as possible.

External fixator systems produced by Ortimplant are produced in different technical specifications, shapes and sizes according to the bone structures and purposes of the treatments performed by medical professionals. These implants and components are compatible with many product structures, and each product structure is designed for specific disease structures. The screws used with these devices do not include cervical vertebra, thoracic vertebra and lumbar vertebra healing products.

External Fixator implants and accessories;

- ✓ Ayfix Dynamic External Fixators & Accessories (1.5/2/2.5/3/3.5/4/4.5/5/5.5/6 mm Schanz screws (or pins)) only for Ayfix External Fixator Device It is used with the set.
- ✓ Mini Rail Fixator & Accessories (1.5/2/2.5/3 mm Schanz screws) are used only with the Mini Rail Fixator Tool Set.
- ✓ Tubular Fixator & Accessories (2.5/3/3.5/4/4.5/5/5.5/6 mm Schanz screws, 5/8/10 mm rods, 5/8/10 mm rod-pin & rod - rod clamps, tubular fixator multiple pin clamps, tubular fixator multiple pin T clamps, tubular fixator oblique/angled supports, tubular fixator 8/10 mm semi-circular rods, tubular fixator compression distraction tube, tubular fixator ankle/elbow hinge) only Tubular Used with the Fixator Instrument Set.
- ✓ Hybrid Fixator & Accessories (4/4.5/5/5.5/6 mm Schanz screws, 155/185/215 mm 2/3 rings, 155/185/215 mm 1/3 rings, hybrid fixator 3-hole slide clamp, washer, hybrid fixator two-hole wire fixation pin, 2.0 mm olive wire, 2.0 mm trocar pointed wire, hybrid fixator independent wire clamp, hybrid fixator additional pin retention clamp, hybrid fixator carbon rod, hybrid fixator wire fixation bolt, hybrid stabilizer straight/angled posts, hybrid stabilizer half pin universal retaining bolt, connecting bolt, nut, studded washer, hybrid stabilizer single pin clamp, hybrid stabilizer hinge clamp, hybrid stabilizer telescopic reduction unit, hybrid stabilizer translational rotation hinge, square socket, cannulated bolt, threaded rod, hybrid fixator extension plate) should only be used with the Hybrid Fixator Instrument Set.
- ✓ Circular External Fixator (Adult) & Accessories (80/100/120/130/140/150/160/170/180/200/220 Half/C/ Foot rings, 90°/120°, 160-180-200 mm femoral belts, multi-pin fixing clamp, oblique support connector, adjustable clamp, telescopic rod, extension rod, double female hinge, male/female half hinge, cannulated/slotted wire fixing bolts, male/female supports, half-pin universal retaining bolt, single fixing bolt, connection bolt M6,

connection/angular/bent/curved plate, square socket, conical/oval/star/slotted washer, M6 nut, threaded/extension slotted rod M6, 1.8/2 mm trocar pointed/bayonet/olive wires , hybrid fixator adaptation for the adult circular fixator) should only be used with the Circular External Fixator (Adult) Instrument Set.

- ✓ Circular External Fixator (Pediatric) & Accessories (70/80/90/100/110/120 mm Half/C/Toe rings, 90°/120° 160-180-200 mm femoral arches, pediatric multi-pin fixation clamp, pediatric oblique support connector, pediatric adjustable clamp, pediatric telescopic rod, pediatric extension rod, pediatric double female hinge, pediatric male/female half hinge, pediatric cannulated/slotted wire retaining bolts, pediatric male/female supports, pediatric half pin universal retaining bolt, pediatric single fixing bolt, connecting bolt M5, pediatric connection/angular/bent/curved plate, pediatric square socket, pediatric tapered/oval/star/grooved washer, pediatric nut M5, pediatric threaded/extension/spline tie rod M5, 1.5/1, 8 mm trocar pointed/bayonet/olive wires, pediatric hybrid fixator adaptation for pediatric circular fixator) should only be used with the Circular External Fixator (Pediatric) Instrument Set.

Specialists are responsible for the consequences of using the wrong implant, not meeting the safety criteria in this user guide, or using inappropriate surgical techniques. In such cases, Ortimplant and its authorized distributor are not responsible to the patient and surgeon.

External fixators and components produced by Ortimplant can only be used with Ortimplant application sets. However, the practicing surgeon should avoid applications using Ortimplant External Fixators and Accessories along with other trademarks of external fixator application kits and sets. Otherwise, the surgeon is responsible for any complications experienced during and after the surgery. Ortimplant and its distributors are not responsible to patients and surgeons in such cases.

During the preparation phase before the surgery, the surgeon must make sure that the external fixator implants, application sets and the implants he will use are compatible with each other. Ortimplant and its distributors are not responsible for the negative consequences of surgical application in case of any complications.

Note: External Fixator implants are manufactured from aluminum, titanium, stainless steel and fiber carbon raw materials..

Note: To achieve the best result; Do not use Ortimplant external fixator implants with parts produced by other companies. Do not reuse external fixator implants due to metal fatigue and risk of infection.

Note to Surgeon: Even if the physician is a knowledgeable intermediary between the patient and the company, the important medical information provided in this document must be disclosed to the patient.

Intended Use

External fixator devices are used for the temporary treatment of orthopedic and traumatological disorders for fixation, distraction and compression of different bones and joints in the human body.

It is mainly used for fixation of the following open fractures:

- ✓ For distraction and compression purposes in open and closed fractures of the femur and in the fracture area,
- ✓ For distraction and compression purposes in open and closed fractures of the tibia and in the fracture area,
- ✓ For open and closed humerus fractures and distraction and compression in the fracture area,
- ✓ Open and closed and/or comminuted fractures of the hand,
- ✓ Open-closed and/or comminuted fractures of the pelvis,
- ✓ For distraction and compression purposes in ulna & radius and fracture area,

It is also intended for the treatment of limb lengthening, deformity and angular correction, malunion, arthrodesis and pseudoarthrosis of long bones.

Indication Information

External fixators and their components are designed to be used in the following treatments;

- ✓ 1. Fractures (especially open and closed fractures),
- ✓ 2. Extremity deformities,
- ✓ 3. Pseudoarthrosis,
- ✓ 4. Joint contractures,
- ✓ 5. Congenital or acquired orthopedic disorders,
- ✓ 6. Arthrosis (arthrodesis),
- ✓ 7. Some joint dislocations,

External fixators and their components can also be used for fixation following deformity correction osteotomies, limb lengthening treatments, replantation treatments, and bone reconstructive procedures.

Note: The size, type and dimensions of the implants and accessories to be used in surgery and the "application sets" suitable for the selected implants and accessories should be determined by the surgeon who will perform the application according to the bone area to be used, the structure of the bone, the condition of the disease, the gender of the patient and the age of the patient. The Ortimeplant catalog contains information about external fixators to be used with implants.

Planning the surgical procedure should be based on detailed evaluation of the patient's x-rays. In this way, the selection of the appropriate external fixator is determined by the surgeon by calculating the size and finish combinations. Pre-operative planning templates are available from the manufacturer. The patient's allergic sensitivity to the implant material should be determined. External fixators and components of different sizes should be available at the time of surgery. Before placing the appropriate implant, the surgeon must ensure that the implant size and type are compared to the label to use the correct size implant. Excess labels inside the box; for reporting and ordering purposes.

Note: External fixator implants, components and application sets suitable for the implant should be selected by the surgeon and her team before the operation, and their operation should be checked and verified.

The surgeon and his team are responsible for selecting the appropriate implant and its components for the patient and the application sets suitable for these implants and components. Ortimeplant and distributors cannot be held responsible for any problems arising from the selection and application of inappropriate implants, components and application sets.

Contraindication Information

The following situations include, but are not limited to, contraindications;

- ✓ Development of active infection or risk of infection (immunodeficiency),
- ✓ Indicators of local inflammation,
- ✓ Increased fever and white blood cell mass,
- ✓ Circulatory disorders and occlusions that prevent adequate blood circulation in the fracture or surgical area,
- ✓ Morbid obesity,
- ✓ Pregnancy (until the patient's suitability for surgery conditions is confirmed by experts),
- ✓ mental illnesses,
- ✓ Some mental or neuromuscular diseases that may cause complications or fixation failure in the postoperative period,
- ✓ Patients with serious congenital anatomical disorders,
- ✓ The presence of another disease or associated use of a medication or device may potentially reduce or prevent the benefit of external fixator implants (e.g., progressive streptococcal infections, congenital functional abnormality, unexplained semimetal elevation, other disabling conditions).
- ✓ Bone stock disorders that prevent fixation or support of the device,
- ✓ Patients with any anomalies in the bone anatomy in the area where the external fixator is intended to be implanted,
- ✓ Rapidly progressive joint disease, bone erosion, osteopenia and/or osteoporosis. Osteoporosis and osteopenia are contraindications that may limit the degree of healing due to decreased stabilization and/or mechanical fixation.
- ✓ Cases of metal allergy to titanium, stainless steel,
- ✓ Situations where bone graft or fusion is unnecessary,
- ✓ In cases where the implant components selected for use are too large or too small for the bone area to be applied,
- ✓ Patients who do not have enough soft tissue to cover the operation area, who do not have sufficient bone stock, or who do not have sufficient bone quality,
- ✓ Patients expected to have low physiological performance,
- ✓ Patients who are not willing to follow postoperative instructions
- ✓ Using the products in situations not specified in the indications,

Note: Due to the high risk of unexpected metal corrosion inside the body, two different metal components should not be used together.

Potential adverse events and complications

All known scientific safety warnings regarding bone surgery will be taken into account by the surgeon and the team during the placement of external fixators and their components. In many cases, negative results may be clinic-related rather than device-related. Potential safety warnings, potential adverse events, and complications related to implantation of external fixators and their components are listed below, but are not limited to;

- ✓ Early or late loosening of the screw components in the external fixator and accessories due to poor bone quality of the patient or incorrect application,
- ✓ Possibility of disassembly, bending and/or breakage of all components within the body.
- ✓ Nail infection may occur when using wire, Schanz pin/screw, Steinman Pin. The patient should be informed by the surgeon. Necessary precautions should be taken with dressing and antibiotic treatment when nail infection occurs.
- ✓ In patients without adequate soft tissue cover, skin penetration due to pressure on the operation side, increase in fibrous tissue, development of skin necrosis and/or pain, and inflammation of surrounding tissues and bone may occur. Improper use of instruments or incorrect placement of the implant may cause damage to surrounding tissue and nerves.
- ✓ After the surgery, changes in bone smoothness, limb weakening and/or length loss in the relevant body part should be monitored by the surgeon and the patient.
- ✓ Development of local and/or general infection,
- ✓ Paralysis (complete or partial), dysesthesia, hypersensitivity, numbness, paresthesia, appearance of radiculopathy and/or pain, depersonalization, loss of neurological function such as neuroma development, loss of consciousness, tingling sensation and/or visible defects,
- ✓ Disorders of bladder fullness, urinary retention and other urinary system disorders,
- ✓ Scar formation can cause nervous system problems or pressure and pain around the nerve,
- ✓ Fracture, microfracture, resorption, development of bone defect or penetration of bone,
- ✓ Pseudoarthrosis, late and poor union of the fractured bone where the surgical operation was performed,
- ✓ Increased/decreased mobility and function in the adjacent joints of the operated bone,
- ✓ The patient cannot continue his daily activities,
- ✓ Loss or decrease in bone density (possibly a result of a decrease in the load on the bone),
- ✓ Complications of the grafted area, including pain, fracture or wound healing problems,
- ✓ Intestinal obstruction or loss of control, gastritis or other gastrointestinal tract disorders
- ✓ deep vein thrombosis,
- ✓ Bleeding, hematoma, seroma, edema, hypertension, vascular occlusion, stroke, excessive bleeding, phlebitis, tissue decay, wound dehiscence, vascular injuries or other cardiovascular system disorders,
- ✓ Avascular necrosis,
- ✓ Shortness in the limb where the implant was applied,
- ✓ Reproductive system disorders, sexual dysfunction, progression of respiratory problems (e.g. pulmonary embolism, pneumonia, bronchitis),
- ✓ Incidence of external fixators into the joint space due to incorrect application,
- ✓ Patients should stay away from strong magnetic fields.
- ✓ The patient should immediately consult his physician if he hears abnormal sounds, pain, or increased temperature in the surgery area.
- ✓ It has been reported that during screwing of the external fixator, especially around the joint, screws of inappropriate size may enter the joints and cause chondral lesions. This should be avoided to prevent chondral lesions and prevent early joint arthrosis. Damage to the epiphysis of long bones during surgery or due to improper placement or incorrectly sized screws and guidewires.

Note: Additional surgical intervention may be required to correct these potential disadvantages.

Intended User Profile

Ortimplant External Fixators should only be used by healthcare professionals.

Surgical procedures should only be performed by medical specialists or medical residents under the responsibility of specialists who have received adequate training in the anatomical structure or formation where the surgery will be performed and are familiar with surgical techniques.

Before performing a surgical procedure, surgical literature on techniques, complications, and hazards should be reviewed. You should carefully read all operating instructions regarding product safety features before use. All necessary materials and documents must be provided by the authorized distributor.

General Safety Warnings

- ✓ Read the User Manual before use and take into account the warnings.
- ✓ The products are NOT STERILE! Before first use, sterilize according to the instructions specified in this manual.
- ✓ Do not use it for purposes other than its intended use.
- ✓ This device is not used to correct any pediatric fractures other than proximal femur fractures.
- ✓ Choosing the appropriate device, appropriate accessory and appropriate application tool set for the patient is very important. Improper use or misalignment of components may cause loss, bending, breakage, or cracking of the device or bones.
- ✓ Until radiographic confirmation of adequate bone development; Other internal or external support devices should be used (orthosis, plaster, etc.) to stabilize the fracture and minimize the load on the fractured area.
- ✓ High pressure occurs in proximal/distal fractures and osteotomies of long bones. This can cause bending or breakage of implant components due to uneven muscle force, which can negatively impact healing.
- ✓ Do not mix incompatible metals.
- ✓ The device is made of metal materials and does not contain waste materials..

Precautions

External fixator implants should not be reused. Loading too early may cause devices to stop working. In order to prevent the implant from causing discomfort during use, check the implants and application sets (in terms of suitability and full functionality) before use. Precautions to be taken before, during and after other procedures are listed below.

Operation

Patient and Implant Selection:

Only patients with indications and no absolute contraindications should undergo surgical procedure. However, the final decision and responsibility belongs to the surgeon who will perform the surgery.

Choosing the appropriate size, shape and design for each patient is important for the procedure to be successful. Metallic surgical implants are subjected to repeated loads during use and are limited in shape, size and design. Unless care is taken in selecting the appropriate implant for the patient, unless the implant is placed properly, and if no precautions are taken to reduce the stress on the implant, such stresses can result in buckling and device loosening, metal fatigue, and eventual fracture before healing. These situations may result in device damage or premature removal.

Before Operation:

- ✓ Only patients who meet the criteria specified in the indications should be selected.
- ✓ The use of external fixators should be avoided in patients who have or are prone to the conditions specified in the contraindications.
- ✓ Care should be taken when transporting and storing external fixators and accessories. Fixators should not be scratched or damaged. They must be protected against aggressive environmental conditions during storage. Implant packages should be checked carefully for damage before use.
- ✓ External fixators of appropriate size and number should be used. Additionally, excess stock and variety should be kept ready during surgery.
- ✓ Surgeons should become familiar with the various components of the external fixator before using the surgical equipment and assemble the external fixator set to verify that all parts and necessary tools are present before the

operation begins. The external fixator and accessories (as described in the description) cannot be completed with components produced by other manufacturers. Different types of metal should never be used together.

- ✓ External Fixators and accessories should be clean and sterile before use. Additional sterile materials (implant size, type alternatives, application sets, hospital orthopedic surgery sets, etc.) should be available to prevent undesirable situations.

During the operation:

- ✓ Extreme care should be taken in the surgical approach to the areas through which the spinal cord, nerve roots and peripheral nervous system pass. Damage to nerves leads to loss of neurological functions. During the surgical procedure, maximum attention should be paid to intracranial, intrathoracic, intrapelvic organs, arterial and venous areas and these structures should be protected. Arterial or venous injuries can result in death.
- ✓ Breakage, slippage, or misuse of instrument sets or implant components may result in injury to the patient or surgical personnel.
- ✓ External fixator implants should not be repeatedly bent excessively. Scratches or nicks on the implant may reduce the functional integrity of the implant; in such cases, due care should be taken by the surgeon to avoid the use of the implant.
- ✓ Whenever possible or necessary, an imaging system should be used to facilitate the surgical procedure.
- ✓ Do not use bone adhesives as it will be difficult to separate external fixators and accessories and their safety and effectiveness have not been established. Heat generated during the healing process can cause neurological injuries and local bone tissue death.
- ✓ Before suturing the soft tissue, change the position of the external fixators if there is an inappropriate external fixator implant position.

After the operation:

The post-operative instructions and warnings given to the patient by the surgeon and the patient complaints associated with them are extremely important;

- ✓ The patient should be informed to avoid placing excessive regional load on the operation area. In such cases, before bone union is completed, the patient should be informed that the external fixator implants may bend and bend, and the fracture position may be distorted. The patient should be informed by the surgeon about the risk of reoperation under these conditions. The surgeon should inform the patient about contraindications that may cause the device to loosen and/or break. If the patient is extremely elderly, frail, or mobile, the risk of bending, loosening, and breaking of the temporary fixation device may be increased during the postoperative rehabilitation period. The patient should be warned to avoid falls and sudden jolts in different positions.
- ✓ To achieve maximum surgical success, neither the patient nor the device should be subjected to mechanical vibration or shock that could cause the device to loosen. The patient should be informed about this possibility and, in particular, his/her participation in any sports activities should be restricted. The patient should be advised to be accompanied by those nearby and to use auxiliary support devices (crutches, walkers, orthoses, etc.) while standing up. The patient should especially be advised to avoid turning movements. The patient should be warned not to drink or consume a product containing nicotine, alcohol, or steroids during the bone healing process.
- ✓ The patient should be warned about the loss of movement in the surgical area and informed about how to compensate for the permanent physical limitation in body movement.
- ✓ Nonunion or delayed union of the fracture; creates excessive and repetitive stress on the implant that stabilizes the fracture. Excessive load may cause the metal device to bend, loosen, or collapse due to fatigue. It is important not to move the surgical area until appropriate bone union is achieved by x-ray examination. If the device is loose, bent and/or warped, it must be replaced or removed without serious injury. The patient should be warned of these dangers and monitored closely until the bone consolidates.
- ✓ Especially in high-risk patients, prophylactic antibiotic use is precautionary before subsequent surgical procedures (such as dental surgery).

- ✓ Implant fixation devices are designed to provide stability at the fracture site during normal fracture healing. Once the bone is fused, these devices have no functional effect and can be removed. The final decision to remove the device is up to the surgeon and the patient. Many patients may experience one or more of the following complications;
 - Corrosion, local tissue reaction or pain,
 - Disruption of implant position (leading to possible injuries),
 - Extra risk of injury due to trauma to the surgical area after surgery,
 - Bending, loosening or breaking of the implant in the surgical area, which may make removal difficult,
 - Pain, discomfort or abnormal movement while the implant is in the body,
 - High risk of possible infection,
 - Bone loss caused by excessive pressure on the bone.

Removing the Implant:

The length of time the device remains on the patient is determined by the surgeon, taking into account the progress of the union and the duration of treatment. The fixator can be removed with a simple outpatient procedure. Before the fixator frame is removed, the patient is given an appropriate dose of paracetamol. The fixing clamps and nuts are loosened and the fixator is separated from the bone screws. Turning the bone screws counterclockwise easily loosens them and pulls them out of the bone. The wounds are closed by placing dry sterile gauze on them.

Note: Ortimplant; recommends removing the implant after the fracture has healed to prevent refracture of the bone and other complications. The external fixator implant and its components should not be reused. It is not possible to use a device removed from one patient on another patient and its use should not even be considered.

Surgical Technique

A separate surgical technique guide is presented to explain the procedure of external fixators during the surgical procedure. The surgeon should be familiar with the surgical procedure before using the products. It is also the surgeon's responsibility to consult with colleagues experienced in the procedures and to become familiar with relevant publications before use. Please review the surgical technique document provided by Ortimplant for the correct application method.

Sterilization

Steam Sterilization Method has been validated and is recommended as the sterilization method for our products. The steam sterilizer must comply with the EN 285 standard and be validated according to the EN ISO 17665 standard.

- ! The steam sterilizer must comply with the EN 285 standard and be validated according to the EN ISO 17665 standard.
- ! It must be ensured that the maximum allowable loading is not exceeded according to the sterilizer manufacturer's instructions.
- ! The properties of the water to be used in the sterilization process must comply with EN 285 and ISO 17665 standards.

Sterilization Method	Temperature	Time
Steam Sterilization	134 °C	5 min

Packaging and Storage Conditions

All Ortimplant External Fixators are cleaned and disinfected according to standards-compliant methods. Cleaning and decontamination involves rinsing with neutral cleaning agents followed by washing with distilled water.

Afterwards, all products are packaged in the clean room.

Note: It should be noted that the packages containing Ortimplant External Fixators are not damaged. If a delivery note or invoice is used, the completeness of all products, including application sets, must be checked before use. Damaged packages and products should not be used and should be returned to Ortimplant or the distributor.

External Fixators should be stored in an environment with humidity below 80% and the environment should not contain corrosive gases.

Disposal Conditions

Disposal of the products must be carried out in accordance with the national legislation of the country where the products are used.

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Manufacturer: Ortimplant Ortopedi ve Tıbbi Aletler Tic. Ltd. Şti.