# CERRAHİ EL ALETLERİ İÇİN GENEL KULLANIM KILAVUZU

## İÇİNDEKİLER

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AAMI TIR 30:2011, bölüm 7.5: A Compendium Of Processes, Materials, Test Methods, And Acceptance Criteria For Cleaning Reusable Medical Devices

ANSI AAMI ST81: Sterilization Of Medical Devices - Information To Be Provided By The Manufacturer For The Processing Of Resterilizable Medical Devices

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EN ISO 17665: Sterilization of health care product- Moist heat-Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices

ISO 15883 serileri: Washer-disinfectors — Part 1: General requirements, terms and definitions and tests

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# INSTRUCTION FOR USE SURGICAL INSTRUMENTS

## CONTENTS

|   |  |
|---|--|
| 1. INTRODUCTION .....   |  |
| 2. PRODUCT INFORMATION .....  |  |
| Intended Purpose .....  |  |
| Intended User Profile .....   |  |
| Product Description .....   |  |
| Warnings and Precautions .....  |  |
| 3. CLEANING .....   |  |
| General Information.....  |  |
| Actions Before Cleaning .....   |  |
| Pre-Cleaning .....  |  |
| Manual Cleaning and Disinfection .....  |  |
| Auto- Cleaning and disinfection by using washer-disinfector.....  |  |
| Drying .....  |  |
| Lubrication/ Assembly.....  |  |
| Inspection.....   |  |
| 4. PACKAGING.....   |  |
| Double Packaging .....  |  |
| Sterilization Containers .....  |  |
| 5. STERILIZATION .....  |  |
| Sterilization Process .....   |  |
| 6. STORAGE.....   |  |
| 7. MAINTENANCE & REPAIR .....   |  |
| 8. COMPLAINTS .....   |  |
| 9. SYMBOLS.....   |  |
| 10. REFERENCES (ONLY FOR SOME COUNTRIES) .....  |  |
| APPENDIX 1- THE CLEANING AGENTS AND DISINFECTIONS USED DURING VALIDATION OF OPERATING INSTRUCTIONS..... |  |
| APPENDIX 2- REPROCESSING OF REUSABLE FLEXIBLE REAMER.....   |  |
| APPENDIX 3- CONTROL GUIDES TO CHECK OPERATE PROPERLY OF SURGICAL INSTRUMENTS.....                       |  |

# INSTRUCTION FOR USE FOR SURGICAL INSTRUMENTS

## 1. INTRODUCTION

- ✚ This document is created according to
  - Medical Devices Regulation (MDR) EU 2017/745
  - ISO 13485 (Medical Devices, Quality management systems, Requirements for regulatory purposes),
  - ISO/TC 170 (Surgical instruments)
  - ISO 7153-1:2016 (Surgical instruments – Materials – Part 1: Metals)
  - ISO 7151:1988 (Surgical instruments — Non-cutting, articulated instruments — General requirements and test methods)
  - ISO 7741:1986 (Instruments for surgery — Scissors and shears — General requirements and test methods)
  - ISO 5832-1:2016 (Implants for surgery — Metallic materials — Part 1: Wrought stainless steel)
  - ISO 12010:2019 (Water quality — Determination of short-chain polychlorinated alkanes (SCCP) in water — Method using gas chromatography-mass spectrometry (GC-MS) and negative-ion chemical ionization (NCI))
  - ISO 13402:1995 (Surgical and dental hand instruments — Determination of resistance against autoclaving, corrosion and thermal exposure)
  - ISO 17664 (Processing of health care products - Information to be provided by the medical device manufacturer for the processing of medical devices,
  - ISO 17665-1:2006 (Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices),
  - ISO 11607-1:2019 Packaging for Terminally Sterilized Medical Devices — Part 1: Requirements for Materials, Sterile Barrier Systems and Packaging Systems,
  - ISO 11607-2:2019 Packaging for Terminally Sterilized Medical Devices — Part 2: Validation Requirements for Forming, Sealing and Assembly Processes, requirements.
- ✚ This document is provided with assembly and disassembly instructions for multi-component surgical handpieces that need to be disassembled before cleaning, and the instructions for use included with the products.

## 2. PRODUCT INFORMATION

### Intended Purpose

- ✚ Surgical instruments are produced to use by surgeons and facilitate to surgeon's work during surgical operation in complex and specific functions such as cutting, dissecting, drilling, pulling, reaming, tearing, clamping, clutching, seeing the surgical site, hemostasia, engraving, scraping, measure the size, probing, draining, aspiration, suturing, ligation, etc.
- ✚ Additionally, to ease implant surgical plates and prostheses, surgical instruments are used.
- ✚ [The surgeon in charge should decide which surgical instruments use to in surgical procedures.](#)
- ✚ It should not be used for any purpose other than surgical procedures.

### Intended User Profile

**Surgical instruments** can only be used by healthcare professionals and only during surgical procedures. Surgical procedures should only be performed by medical specialists who have had sufficient education in anatomical structure or formation in which the surgery will be performed and who are familiar with surgical techniques, or medical residency students under the responsibility of expert. In addition, veterinarians may also use some surgical hand tools (the indication of which will be determined by themselves) in surgical procedures belonging to their own branch. Before performing a surgical procedure, should be reviewed to surgical literature about techniques, complications, and hazards; should be read carefully all instructions for use about product safety specifications before using. All materials and documents required should be provided by the authorized distributor.

### Product Description

#### Material Used for the Product

**Surgical Instruments,**

✚ **ISO 7153-1:2016** (Surgical instruments – Materials – Part 1: Metals) It is produced from medical grade stainless - heat resistant steel, titanium and aluminum and alloys of these metals in accordance with the standard.

✚ All surgical instruments; Manufactured from **ISO 5832-1** stainless steel. The materials used for each product can be seen in the product catalog and product list

✚ Stainless steel; Due to its chemical (high corrosion resistance and biocompatibility) and physical (tensile strengths; stainless steel's physical properties of fractures are the best) properties, they are mostly used in the manufacture of surgical instruments and are validated raw materials.

✚ Some parts of the instrument, such as the handle or body, are made of **aluminum, fiber, silicone or plastic** for better holding of the instrument, applying force, and taking measurements during the surgical procedure.

#### **(Table 1)**

✚ Steel; It is fully biocompatible, corrosion resistant and non-toxic to the human body.

✚ **Reusable surgical instrument groups** are shown in **Table 2**. Great attention must be paid to cleaning, disinfection and sterilization procedures for reusable instruments.

✚ **Disposable surgical handpiece groups** are shown in **Table 3**. Since this group of hand tools are not designed to work as intended after first use; should not be reused. Reuse may result in reduced reliability, performance and/or compatibility of the instrument. In this case, the responsibility belongs to the user (responsible surgeon).

✚ **Surgical instruments** are sold non-sterile.

✚ However, it is in critical class according to ISO/FDA guidelines and Spaulding Classification Chart and should be used sterile during procedures. For this reason, before the first use and every subsequent use, **ISO 17664:2017** (Processing of health care products - Information to be provided by the medical device manufacturer for the processing of medical devices) and **ISO 17665-1:2006** (Sterilization of health care products) - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices) and should be cleaned, disinfected and sterilized according to the procedures specified in this document.

✚ Surgical instruments do not implant.

Use of surgical instruments other than the above-mentioned indications may result in damage or breakage of the instruments. Proper cleaning, handling, sterilization, and standard routine maintenance (such as sharpening, sharpening) will ensure that operate surgical instruments as intended and extend their lifetime.

Inspect tools for possible damage, wear, corrosion, discoloration, or non-functioning parts before use. Carefully inspect critical, inaccessible areas, connections, and all moving parts. Damaged, incomplete, and inoperable surgical instruments should not be used. Contact your local sales representative for repair or replacement.

## **Surgical Instruments Product Groups**

### **1. SIs with Plastic, Polyethylene, etc. parts**

1. Screwdrivers
2. Retractors
3. Drill Guides

### **2. Disposable SIs**

1. Drill bits (4.5 mm and thinner ones)
2. Gigli saw wire

### **3. Class I/II SI Groups According to Intended Purpose**

1. Cutter (Dissection) SIs
  - a) Scalpel Handles (Lancets)
    - a. Operating scalpel handles
  - b) Blades
    - a. Fistula knife
    - b. Tonsil knife
    - c. Resection knife
  - c) Scissors
    - a. Bandage scissors

- b. Dissecting scissors
- c. Iris scissors
- d. Mayo scissors
- e. Metzenbaum (Lahey) scissors
- f. Micro scissors
- g. Nail scissors
- h. Operating scissors
- i. Suture scissors
- j. Vascular scissors
- d) Trocars
- e) Rongeurs
- f) Saws
  - a. Bone saws
    - i. Powered saws
    - ii. Unpowered saws
- g) Diathermy/Cautery
- h) Elevators
- 2. Holder (Grasper, Clamper and Occluding) SIs
  - a) Clamps
    - a. Bulldog clamps
  - b) Hemostatic clamps
    - a. Mosquito clamps
  - c) Dissection clamps
  - d) Forceps
    - a. Tissue forceps
    - b. Dissecting forceps
    - c. Ligature forceps
    - d. Artery forceps
    - e. Mosquito forceps
    - f. Hemostatic forceps
    - g. Bronchoscopy forceps
    - h. Biopsy forceps
    - i. Mixer forceps
    - j. Smoothed forceps
    - k. Toothed forceps (to grasp skin)
    - l. Allis forceps (for colporrhaphy (anterior and posterior repair))
    - m. Babcock forceps (to grasp intestine and towels)
    - n. Stone forceps
  - e) Needle holders
    - a. Needle holder
    - b. Microvascular needle holders
- 3. Retractor (Dilatator, Exposing) SIs
  - a) Automatic (self-retaining) retractors
    - a. Adson self-retaining retractor
    - b. Balfour abdominal retractor
  - b) Retractors with blade
    - a. Retractor blade
  - c) Retractors without blade
    - a. Wolkman retractors
    - b. Kocher retractors
    - c. Langenbeck retractors
    - d. Army navy / Farabeuf retractor
  - d) Hooks
    - a. Skin hooks

- b. Bone hooks
  - c. Nerve/vessel hooks
  - d. Tendon hooks
- e) Spatulas
- f) Speculums
  - a. Vaginal speculum
- 4. Reamer, Driller SIs
  - a) AWL, drilling tip
  - b) Reaming Scrapper
  - c) Trocars
- 5. Scrapper SI
  - a) Curettes
  - b) Scrapper
  - c) Cutter and gujs
  - d) Spoons
- 6. Probes & Scopes
  - a) Double Ended Probe
  - b) Urethral sounds
  - c) Laryngoscopes
  - d) Anoscopes
  - e) Proctoscopes
  - f) Sigmoidoscopes
- 7. Suction & Aspiration SIs
  - a) Suction tubes
  - b) Tracheal tubes
- 8. Measurement SIs
  - a) Tape measure of oil cloth
  - b) Rulers
  - c) Calipers
- 9. Hammers SIs
  - a) Hammers
- 10. Markers SIs
- 11. Other SIs
  - a) Tromels
  - b) Instrument trays
  - c) Bathtub cover
  - d) Cotton box
  - e) Laboratory dish
  - f) Kidney tray

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## References

1. [https://en.wikipedia.org/wiki/Instruments\\_used\\_in\\_general\\_surgery](https://en.wikipedia.org/wiki/Instruments_used_in_general_surgery)
2. <http://research.sklarcorp.com/4-basic-types-of-surgical-instruments>
3. Heller, Michelle. (2016). Clinical Medical Assisting: A professional, Field Smart Approach to the Workplace.
4. 2 Surgical Instruments slide share
5. <https://www.surgicalholdings.co.uk/history-of-surgical-instruments.html>

## Warnings and Precautions

- ⚠ A disposable medical device; when provided nonsterile and if the customer need sterilize before use, the relevant section of these instructions should be performed as long as have not any specific instructions in the prospectus.
- ⚠ By reason of the fact that they are not designed to function as intended after the first use, **do not re-use disposable surgical instruments**. Re-use may result in a decrease in the reliability, performance, and/or compatibility of the tool. In this case, the responsibility belongs to the user(responsible surgeon). Check the label of the device to specify one or more than one and/or cleaning and re-sterilization process.
- ⚠ Visually inspect surgical instruments whether any damage or abrasion. Instruments with broken, fractured, or eroded parts/surfaces should not be used, should be replaced immediately.
- ⚠ Check the cutting edges of the drill bits are smooth and continuous, it does not contain large cracks or chips that adversely affect cutting performance.
- ⚠ Verify that the matching surfaces are working purposefully and without complications. Do not use surgical instruments that do not perform satisfactorily.
- ⚠ Be careful when using sharp and pointed tools to prevent injury.
- ⚠ Do not normally use in surgical procedures when in touch with infective texture including low or high risks *Contagious Spongiform Encephalopathies* (e.g.: *patient diagnosed with or suspected of Creutzfeldt Jakob Disease (CJD)*). In that case, does not apply procedures of sterilization. The product should be treated as medical waste.

## 3. CLEANINGS

### General Information

- ⚠ In this section, Information on 2 methods used to clean surgical instruments is presented.
  - Manual Method
  - Method Using an Automatic washer-disinfector
- ⚠ Complete manuel and automatical cleaning totally before the sterilization. Automatic methods may not clean enough the instruments alone.
- ⚠ However, the automatic wash method is more practical and safe. Also, the person who is working on this process is less exposed to contaminated instruments and used cleaning agents. Therefore, the automatic method should preferably be used after removing coarse dirt and residues with manual cleaning.
- ⚠ The cleaning solutions should be changed often and without being allowed to become very contaminated.

- ✚ In the cleaning process (whichever method is used); The employee must use appropriate personal protective equipment (PPE) and equipment at all times.
- ✚ Pay particular attention to the instructions on the correct handling and use of the product provided by the cleaning agent manufacturer.
- ✚ Only specially formulated cleaning agents and/or disinfectants (detergents) should be used in the cleaning and disinfection of surgical instruments. Guidelines for immersion in cleaning solutions and/or disinfectants given by detergent manufacturers must be followed.
- ✚ Selection criteria of the relevant detergents are presented in **ANNEX-1**.
- ✚ The quality of the water used in cleaning is important. Freshly prepared pure water or ultrapure water application is recommended. Mineral residues in hard water may cause staining of the surgical instrument. In addition, it can prevent effective cleaning and decontamination with a high risk of contamination containing microorganisms and endotoxins.
- ✚ Special instructions for flexible reamers are given in **ANNEX-2**.

Trays are for the sterilization, transportation, and storage of surgical instruments. These trays are not designed for cleaning and disinfection in a fully equipped state. Instruments must be removed from the tray to obtain satisfactory cleaning results.

### Jobs Before Cleanings

- ✚ Bring surgical instruments to the place which is doing cleanings as soon as possible. If the delay possibility of transfer to transact place is high, cover surgical instruments with a damp cloth or keep the instruments in closed boxes to avoid the dirt from drying out.
- ✚ Remove dirt with a paper towel shortly after the surgical process (within max 2 hours). Also, wash surgical instruments liberally under warm to cold (43 ° C / 110 ° F) tap water or transfer to a bath with an aldehyde-free disinfectant solution that meets the criteria given in Annex-1. Do not use warm water, alcohol, disinfections, antiseptics, saline, or chlorinated water.
- ✚ Take care to avoid mechanical damage to surgical instruments during transport. For example, clean heavy instruments, precision instruments, and cutting/drilling instruments in separate baskets.
- ✚ Similarly, do not apply the same cleaning action to instruments of different metallic compositions. For example, the cleaning procedures are not the same for surgical instruments made of plastic, stainless and heat-resistant steel, titanium, and aluminum. To protect the surface coating of ebonized instruments, keep them separate from other instruments and use a more delicate cleaning procedure (avoid mechanical cleaning and abrasive cleaners).
- ✚ Colored anodized aluminum tools can lose their color in conventional and mechanical cleaning processes.
- ✚ Disassemble the device if needed for cleaning. For this, see or request the instruction manual and the instructions in the surgical technique or additional information available.

Pay attention to the sharp and sharp edges of surgical instruments during surgical procedures, as well as during handling, cleaning, disinfection, sterilization, and packaging, to prevent worker injury and instrument damage. Put your health and safety first and protect yourself.

Microsurgery, plates, and sensitive instruments should be cleaned chemically or manually. Should not be put into the ultrasonic cleaner. In addition, care should be taken to protect sensitive microsurgical instrument tips during the entire cleaning and sterilization process.

### Pre-Cleanings

**Warnings:** In cases where direct cleaning and disinfection are involved, the pre-cleaning step can be skipped. It is recommended to pre-clean highly contaminated surgical instruments in an ultrasonic bath for automatic cleaning.

### Essential Equipment:

- ✚ Cleaning bath/case which is large enough accommodate all instruments,
- ✚ A fresh cleaning solution containing a cleanser for manual cleaning ( temperature should not exceed 50 ° C) that matches the criteria given in **Annex-1**,
- ✚ Cleaning wires for smooth/hard brushes, bottle brushes, and cannulations, etc.
- ✚ Personal protective equipment (at least work uniform, glove, face/eye shield)
- ✚ Absorbing paper
- ✚ Injector (1-50 ml volume based on the magnitude of the cannulation to be wash)

### Procedure

- ✚ Immerse surgical instruments in the solution of the cleanser. Remove accumulated dirt with a paper towel and cleaner solution. Open the mouthpieces of hinged instruments.
- ✚ Make sure that reaching the solution to all surfaces of the instruments, cannulation, etc., and no air is trapped inside the instruments. If necessary, use an injector for this.
- ✚ Clean thoroughly surgical instruments using with appropriate brushes or cleaning fibers (generally only use soft brushes, do not use metal brushes or steel wool)
- ✚ For bone-cutting features such as drill bits, reamer threads, and drill teeth, use a firm bristle brush.
- ✚ Use a suitable diameter bottle brush for cannulations. Make sure the brush passes the full length of each cannulation at least three times.
- ✚ Rinse under running water for at least 1 minute until all traces of cleaning solution are removed.
- ✚ Pay particular attention to cannulations and blind holes, hinges, and connections between mating parts. These should be thoroughly rinsed at least three times with a syringe (volume 1-50ml).
- ✚ Visually inspect for dirt and repeat the above steps if necessary.
- ✚ Wait for the instruments to float on absorbent paper or proceed to the cleaning step immediately.

## **Manual Cleaning and Disinfection**

### **Cleaning**

#### **Required Equipment:**

- ✚ Ultrasonic bath large enough to completely immerse the surgical device (frequency 25-50kHz, temperature - according to detergent manufacturer's instructions)
- ✚ Detergent (low foaming and protein solvent etc.) and rinse aid designed for manual cleaning and suitable for ultrasonic cleaning that meets the criteria given in **Annex 1**.
- ✚ Suitable brushes (only soft brushes, never metal brushes or steel wool) or cleaning wires (for small channels) to reach all parts of the device.
- ✚ Syringe (in a volume of 1 - 50 ml depending on the size of the channels to be washed)
- ✚ Freshly prepared pure water, ultrapure water, or sterile water for rinsing purposes

#### **Procedure:**

- ✚ Prepare an ultrasonic bath containing cleaning solution at the concentration and temperature specified in the detergent manufacturer's instructions.
- ✚ Immerse the device completely and keep the bath active for at least 15 minutes.
- ✚ Thoroughly clean the device using suitable brushes or cleaning wires (for small channels) and paying attention to rough surfaces and lines where dirt may not come off with the brush. In addition, pay special attention to cannulations and blind holes, joints between hinges, and mating parts. Rinse cannulations at least three times with a syringe.
- ✚ After the cleaning step in the ultrasonic bath is completed, if there is any encrusted dirt on the device that needs to be brushed off, repeat the cleaning step as above.

### **Disinfection:**

#### **Required Equipment:**

- ✚ Large bath, whose temperature is under the detergent manufacturer's instructions, large enough to completely immerse the surgical device.
- ✚ A disinfectant designed for manual disinfection that meets the criteria given in **Annex 1** and compatible with the applied cleaning detergent.
- ✚ Syringe (in a volume of 1 - 50 ml, depending on the size of the channels to be flushed).
- ✚ Freshly prepared purified water, ultrapure water, or sterile water for rinsing purposes.
- ✚ Filtered medical compressed air (if available) or clean, lint-free disposable wipes.

#### **Procedures:**

- ✚ Prepare a bath with a disinfectant solution at the concentration and temperature specified in the detergent manufacturer's instructions.
- ✚ Immerse the appliance completely, at least for the time specified in the detergent manufacturer's instructions.
- ✚ Rinse under the running water of the specified quality for at least 1 minute until all traces of disinfectant solution are removed. Pay special attention to cannulations and blind holes, hinge, and connections between mating parts. Rinse at least 5 times with a syringe (volume 1-50 ml).
- ✚ Dry the surgical device with medical compressed air, clean, lint-free, disposable wipes (if necessary, allow to dry for a maximum of 2 hours in a clean area) or by heating in an oven below 110 ° C.

- ✚ Visually inspect and, if necessary, repeat the entire manual cleaning and disinfection process.

The decontamination/disinfection process does not sterilize surgical instruments. Remember that it needs to be treated separately as outlined in the STERILIZATION section.

### **Automatic cleaning and disinfection using a washer-disinfector (recommended)**

#### Required Equipment:

- ✚ Washer-disinfector properly assembled, proper, regularly maintained, and tested, with the required approved efficacy (e.g. CE mark or FDA approval under ISO 15883 series)
- ✚ Filtered air for adequate rinsing steps and active drying program (application of rinse aid is not recommended, the danger of residue) with approved thermal disinfection program (A0 value > 3000 - or for older hand tools - at least 5 minutes application at 90 ° C; chemical disinfection program is not recommended - the risk of disinfectant residue on the device)
- ✚ Final rinse/disinfection with freshly prepared distilled water or ultrapure water
- ✚ Detergent and rinse aid designed for automatic cleaning that meets the criteria given in **Annex-1**.

#### Procedure:

- ✚ For automatic washer-sterilizers, always follow the manufacturer's specifications and use a free-rinse, low-foaming detergent with a neutral pH (6.0 - 8.5).
- ✚ Put surgical instruments in the washer-disinfector.
- ✚ Connect the cannulations to the rinse ports of the washer-disinfector. If direct connection is not possible, position the cannulations directly on the injector jets or on the injector collars of the injector basket.
- ✚ Prevent instruments from touching each other (Movement during washing may scratch and damage anodized instruments or the washing process may stop)
- ✚ Arrange surgical instruments so that cannulations are not horizontal and blind holes slope downward
- ✚ Connected instruments such as scissors must be in an open position.
- ✚ Run the washer-disinfector program (loop)
- ✚ When the cycle is complete, drain the washer-disinfector
- ✚ Visually inspect each tool for dirt residue and dryness.
- ✚ If dirt remains, repeat the cleaning process along with the pre-cleaning phase. The remaining wetness can be removed with medical-grade compressed air, clean, lint-free, disposable wipes (leave to dry for a maximum of 2 more hours in a clean area if necessary), or by heating in an oven below 110 ° C.

Chemical disinfection programs are not recommended due to the possibility of chemical residues remaining on the instrument. These residues may interfere with the effectiveness of sterilization.

### **Drying**

- ✚ Before instruments are packaged and sterilized, they must be thoroughly dried and all remaining moisture removed. Use a soft, absorbent towel/cloth to dry the exterior surfaces. Compressed air or 70% alcohol rinsing can additionally be used to aid the drying process.

### **Oiling/ Assembly**

- ✚ Lubrication is important for moving surgical instruments that consist of more than two parts. Particular attention should be paid to the lubrication of joints, locks, and moving parts. Oil dry instruments only.
- ✚ Reassemble instruments as needed before attaching them to baskets or trays. Examine instruments for bent tips, dents, cracks, misalignment, and corrosion. Remove stained, discolored, or damaged instruments. Mechanically test the working parts to verify each tool is working correctly.

Do not use mineral oil, petroleum, or silicone-based products. Use a silicone-free, water-soluble lubricant before sterilization to lubricate the box locks and connections.

### **Inspection**

Should be checked all surgical instruments before sterilization preparation. For this, visual inspection without using a magnifying glass under sufficient lighting conditions will be sufficient. All parts of instruments should be inspected for visible dirt and/or corrosion.

- ✚ Visually inspect devices for damage or wear. Instruments with broken, cracked, chipped, or worn parts or surfaces should not be used and should be replaced with new ones.
- ✚ Check that the tips and sharp edges of the drills are smooth, uninterrupted, and free of large cracks or chips that could adversely affect cutting performance.
- ✚ Verify that the mating surfaces are functioning as intended and that the powering device interfaces are uncomplicated.

- ✚ In addition, the following observations should be paid attention to:

- Places that collect dirt such as joints, hinges, shafts of flexible reamers,
- Recessed features (holes, cannulations, etc.),
- Features where dirt can be pressed in contact with the device, such as drill flutes adjacent to the cutting tip, sides of the teeth of the drills and rasps,
- Cutting edges should be checked for sharpness and damage.

✚ [In addition, the following functional checks should be made when appropriate.](#)

- ✚ Mountable instruments should be checked for proper assembly.
- ✚ Surgical instruments with moving parts should be practiced to check correct operation.
- ✚ Rotational instruments (drill bits, reamers, etc.) should be checked for straightness. (check instruments by simply rotating them on a flat surface).
- ✚ Flexible instruments (intramedullary reamers, etc.) should be checked for damage.

Note: The cycle (continuity) of surgical instruments' reuse and sterilization has an impact on their lifespan. The period of use of surgical instruments depends on many factors, including the healthcare professional using it, the methods and duration of each use, and the handling/handling between this use. Therefore, it does not specify the maximum number of uses normally suitable for reusable surgical instruments. The best way to determine whether the surgical instrument is at the end of its useful life is to carefully examine the instrument and perform a functional test before use.

✚ [See Annex 2 for more details.](#)

**Figure 1. (A)** A bone/nail cutter with a broken mouth and **(B)** a fractured and corroded osteotome

## 4. PACKAGING BEFORE STERILIZATION

### Double Packaging

✚ Surgical instruments; if appropriate, should be mounted on special trays. Trays should be double wrapped. Drying time can be extended up to 45 minutes for trays containing large polymer components. Finally, the following requirements must be met for the packaging of sterilized surgical instruments:

- EN ISO 11607 standard series,
- Suitable for steam sterilization
- Adequate protection of instruments and sterilization packages against mechanical damage,
- Grading according to the weight of the instrument and implant tray.

### Sterilization Container

✚ In addition to the commonly used double-pack sterilization method, it recommends the use of containers

for the sterilization of reusable surgical instruments.

**Note:** Depending on loading and packaging configurations, a modified drying process may be required for detachable chamber trays. It recommends increasing the drying time to 45 minutes or separating the detachable chamber trays into two parts and placing them side by side in full-size containers. It is the responsibility of the hospital to verify the respective drying time with the sterilization equipment used. For the use of Sterilization Containers, please see Sterilization Container Systems Instructions for Use.

## 5. STERILIZATION

✚ Sterilization of instruments can be carried out with steam (autoclaves) or ethylene oxide (EtO) gas. However, it recommends steam sterilization processes. Surgical Instruments are sterile after steam sterilization (autoclave); however, it should be borne in mind that the design and performance of the autoclave can affect the effectiveness of sterilization. Various indicators etc. in each sterilization process it is necessary to verify that the instruments are sterile.

Surgical instruments containing thermolabile materials should not be subjected to additional loads in the autoclave.

### Sterilization Process

✚ It recommends steam autoclave sterilization with pre-vacuum (air extraction from the autoclave) function. Autoclaves must be validated/calibrated, maintained, checked, and comply with their requirements under EN 285 / EN 13060, EN ISO 17665, and ANSI AAMI ST79.

✚ Surgical instruments should be packaged in an open-mouthed position or disassembled, if appropriate. The instruments can be packaged in rigid containers, woven or non-woven materials. The packaging must ensure the sterility of instruments and allow removal without contamination until the packaging is opened for use in the sterile field.

✚ An autoclave cycle for the sterilization of complete surgical instrument containers/trays. Instruments should be sterilized in the mounting condition is stored in the tray. That is, if the brackets or recesses in the tray are designed to accommodate multi-component instruments in their assembled state, these instruments do not need to be disassembled for sterilization. [The process parameters shown below are recommended by](#)

✚ The ultimate responsibility for verifying sterility using the parameters provided by the equipment and processes of the healthcare facility rests with the healthcare facility. To ensure optimum processing, all cycles and methods should be validated for different sterilization chambers, packaging methods, and/or various loading configurations.

✚ [An exemplary sterilization processes](#)

| Sterilization Method                                | Temperature             | Minimum Exposure Time   |   |
|---|-------------------------|---|---|
|   |                         | Packaging (wrapped)   | Container System  |
| Prevacuum   | 132°-135°C (270°-275°F) | 4 mins  | 4mins (classic or bottom perforated)                                    |
| Gravity   | 121°-123°C (250°-254°F) | 15-30d mins   | 40 mins (only bottom perforated)  |
|   | 132°-135°C (270°-275°F) | 10-25 mins  | 30 mins (only bottom perforated)  |
| IUSS Pre-vacuum (Immediate use steam sterilization) | 132°C/270°F             |   | 3 mins (non-porous products)<br>4 mins (non-porous and porous products) |
| IUSS Gravity  | 132°C/270°F             |   | See manufacturer's instructions for use                                 |
| Ethylene Oxide (EtO)                                | 55°C/131°F              | Concentration: 725mg/L<br>Exposure Time: 60 minutes<br>Humidity: ≥ 50%<br>Aeration: 8 hours (minimum) |   |

✚ An alternative sterilization process

|               |  |
|---------------|--|
| Method        | Steam heat sterilization under EN ISO 17665 or ANSI/AAMI ST79  |
| Loop          | Pre-Vacuum or steam saturated with fractional compressed air   |
| Exposure Time | 3-4 mins (WHO and Robert Koch Institute etc. Can be extended up to 18 minutes to comply with the recommendations.) |

|                       |  |
|-----------------------|--|
| <b>Temperature</b>    | 132-137 °C (270-277 °F) or 134-138 °C (273-280 °F) |
| <b>Drying Time</b>    | 30 min. (minimum, reservoir)                       |
| <b>Cooling Period</b> | 60 min. (minimum, at ambient temperature)          |

**Exposure time:** The time that the load and the entire room are at sterilization temperature

**Drying time:** The time during which the steam is removed from the chamber and the room pressure is reduced to allow evaporation of the condensate from the load by prolonged evacuation or injection and extraction of hot air or other gases. Since the drying time varies depending on the configuration of the load, the method of wrapping, and the material, the healthcare facility must verify the appropriate drying time with the sterilization equipment used.

- For successful sterilization, the autoclave , sterilization containers, and packaged sterile material (such as surgical instruments and implants) must be free of water/moisture residues. A wet package and excessive moisture can create convenient environments for the growth of microorganisms and cause biological contamination problems.

- Therefore, before the surgery in the operating room, it should be checked whether there is water/moisture at the bottom of the container, and whether the surgical cloths are moist, together with the control of biological and chemical indicators inside the sterilization container.

- Autoclave problems (for example, errors in the pressure and vacuum system, clogged drain filter, etc.) and/or operator errors (closing the filter part of the container with material, filling the autoclave with too much material, closing the discharge part of the autoclave, choosing the wrong sterilization cycle, etc.) The manufacturer and its authorized representatives are not responsible for surgical hand tools that are not sterile, corroded and/or stained as a result.

Disposable surgical instruments such as drill bits and Gigli saws should not be reused. It does not recommend the use of "high speed" sterilization for reusable instruments.

## 6. STORAGE

✚ After sterilization, surgical hand tools in sterilization packages should be stored in a dry and dust-free environment.

✚ Shelf life of surgical instruments; Depending on the sterile barrier, storage type, environmental, and usage conditions, it can be extended up to 1 year.

✚ The maximum shelf life for sterilized surgical instruments should be defined by each healthcare facility before use.

✚ Storage conditions in the warehouse should be explained. Look at the information on the label. It should be stored at room temperature and humidity should not exceed 80%.

✚ Surgical instruments are manufactured from medical grade stainless - heat resistant steel, titanium and aluminum and alloys of these metals. Due to the stainless steel content, the device should be kept away from high magnetic fields.

✚ In addition, the device can conduct electricity and therefore should be kept away from external electrical contacts.

✚ Surgical instruments should not be inadvertently interacted with electrocautery. There is a risk of cautery burn as a result of electrical conduction.

✚ The devices are not affected by radiation, they are radiopaque, they do not accumulate radiation.

✚ Hand tools should be stored at room temperature and at a humidity not exceeding 80%.

✚ Since they are conductive, the operating room should not be used before it is completely cooled.

## 7. MAINTENANCE AND REPAIR

✚ Moving surgical instruments, such as scissors and clamps, which consist of multiple parts, may experience corrosion and wear over time due to friction. In order to reduce this problem, hand tools should be regularly maintained and repaired when necessary.

✚ If your Surgical Instruments require repair or maintenance, you can use sponge, foam, bubble wrap, etc. to protect the instruments. wrap it with a protector and return it in a sturdy box. Send packaged instruments to:

✚ For instruments returned to repair; user declaration stating that each instrument has been thoroughly cleaned and disinfected. Failure to provide proof of cleaning and disinfection will also result in a cleaning fee and delay in the repair of your appliance.

✚ Products repaired are guaranteed for 90 days for workmanship and parts when normally used for surgical purposes. Any workmanship or part proven to be defective will be replaced or repaired at authorized discretion, free of charge to the customer.

## 8. COMPLAINTS

✚ The healthcare professional who has any complaints about surgical hand tools or is not satisfied with the quality, identity, durability, reliability, safety, functionality, and/or performance of the product should inform the authorized.

✚ In the event of failure and/or damage to the tool or serious injury or death, the authorized should immediately in writing, by telephone, fax, etc. must be informed with. When filling the complaint form, please fill in the reason for the complaint section. Indicate the serial number and lot number of the instrument (s), a description of the complaint, your name, address, and whether you would like a written reply report for the distributor.

### More Information:

✚ For additional info about this system or any other complaint, please get in contact that is mentioned communication address above.

## 9. SYMBOLS

Lot Numbers

Model no / reference no (Catalog number)

Production Time

Product Barcode

Consult the user manual for warnings and precautions that cannot be presented on the label.

Humidity limitation symbol: Store in environments not exceeding 80% humidity

Do not re-use.

Do not use it if the package is damaged.

Read the user manual

Non-sterile

## 10. REFERENCES (FOR SOME COUNTRIES ONLY)

AAMI TIR 30:2011, Chapter 7.5: A Compendium of Processes, Materials, Test Methods, And Acceptance Criteria for Cleaning Reusable Medical Devices

ANSI AAMI ST81: Sterilization of Medical Devices - Information to Be Provided by The Manufacturer for The Processing of Re-sterilizable Medical Devices

ANSI/AAMI ST77: Containment devices for reusable medical device sterilization

ANSI/AAMI ST79: Comprehensive guide to steam sterilization and sterility assurance in health care facilities

EN 10270-3 / EN 10088-1: Steel wire for mechanical springs. Stainless spring steel wire

EN ISO 11607: Packaging for terminally sterilized medical devices

EN ISO 17664: Processing of healthcare products- Information to be provided by the medical device manufacturer for the processing of medical devices

EN ISO 17665: Sterilization of health care product- Moist Heat-Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical devices

ISO 15883 series: Washer-disinfectors — Part 1: General requirements, terms and definitions, and tests

World Health Organization (WHO) and Robert Koch Institute (RKI) recommendations on exposure time: WHO / CDS / CSR / APH / 2000.3, Annex III, part 2, 6. RKI "Bundesgesundheitsblatt" Communication, April 2002, 45: 376-394, table 11

## **ANNEX-1**

### **CLEANING AGENTS AND DISINFECTANTS USED DURING CLEANING VALIDATION**

✚ pH neutral enzymatic detergents should be used for cleaning and/or disinfection of surgical hand tools made of aluminum alloy. Contact with strong alkaline detergents, iodine, or chlorine-containing solutions should be avoided as they may react chemically with aluminum and cause damage to the handpiece.

✚ For automatic cleaning (washer-disinfector), medium alkaline detergents and brighteners (up to pH 10.9) should be preferred. However, it is recommended to use only the substances recommended for use in these machines. In all cases:

- Follow the indications, instructions, and warnings provided by the cleaning agent and/or disinfectant supplier,
- Choose only detergents designed for cleaning and/or disinfection of surgical instruments made of metal and plastic, and
- Choose only disinfectants with approved efficacy (look for approval or marking such as CE, FDA, VAH/DGHM, etc.).

**Make sure that the items listed below are not the material of the selected cleaning or disinfection detergent:**

- Organic, mineral, and oxidizing acids (minimum accepted pH value 5.5)
- Strong lye (maximum accepted pH value 10.9)
- Organic solvents (acetone, ether, alcohol, gasoline, etc.)
- Oxidizing agents (peroxides, hypochlorite, etc.)
- Halogens (chloride, iodine, bromine, etc.)
- Flavored, halogenated hydrocarbons

## **ANNEX-2**

### **RECONSTRUCTION OF REUSABLE FLEXIBLE REAMERS**

#### **Reprocessing:**

Reusable flexible reagents are considered to be labor-intensive hand tools in terms of cleaning. In any case, cleaning, sterilization, inspection, and maintenance instructions must be followed strictly and precisely. In addition, the following technical advice should be included in practice:

✚ The accumulated dirt should be removed with paper towels immediately after application. In addition, it is recommended that handpieces be washed or transferred to a bathtub containing an aldehyde-free disinfectant solution. If the dirt is not removed immediately after the shafts are used, there is a possibility that incrustation will occur, which will permanently swell in subsequent sterilization.

✚ The flexible reamer shaft needs to be flexed and rotated during each brushing or washing to facilitate dirt removal.

✚ The automatic cleaning process should be preferred to the manual method and should be used as much as possible. The automatic cleaning process has higher manufacturability and reliability.

✚ **Automatic method:** Connect the cannulation to the rinse ports of the washer-disinfector. If direct connection is not possible, position the cannulations directly on the injector jets or on the injector collars of the injector basket.

✚ When the manual method is applied, the pre-cleaning step must be fulfilled. It is suggested that this step be carried out in an ultrasonic bathtub.

**Preventive Maintenance:** The best way to determine whether the surgical instrument has reached the end of its useful life is to carefully examine the instrument and perform functional tests before preparing for sterilization.

In particular, the spiral element should be visually inspected for damage, as the deformed spiral will create a risk for the cleaning of the tool and cause dirt to enter the spiral.

✚ Flexible reamer shafts are manufactured using fixed stainless device steels in a controlled process environment and are welded to the flexible reamer shafts under inert gas conditions to prevent oxidation (corrosion) throughout their lifetime.

#### **Description of Cleaning and Disinfection Process:**

Flexible reamer shafts were subjected to standardized microbiological efficacy control studies to characterize the described automatic and manual methods.

Shafts contaminated with spore-containing artificial dirt in worst-case conditions. After the shafts were reworked according to the instructions, the reduction in spores was determined. Manual and automatic methods have been described separately. A reduction of at least 3 log steps in the number of spores was

considered acceptable.

### **ANNEX-3**

#### **GUIDES TO CHECK THAT SURGICAL INSTRUMENTS WORK PROPERLY**

The following guidelines should be applied to all surgical instruments labeled for multiple uses. All functional checks and inspections described below also cover their interfaces with other instruments or components.

The failure modes below may be caused by the end of life of the product, misuse, or improper maintenance.

It does not specify the maximum number of uses for normally reusable Surgical instruments. The useful life of these handpieces depends on many factors, including the method and duration of each use and the handling between that use. The best way to determine if the surgical device is at the end of its useful life is to carefully examine the device and perform a functional test before use. However, for certain instruments, end-of-life is defined, validated, and specified either by the number of uses or by the expiration date.

#### **Functional check for multi-use Drill Bits**

Description and Function: Multiple use drill bits, cannulated drill bits, drills, plugs, core drills Possible Failure Modes:

- ✚ Defective coupling end (worn)
- ✚ Dull and non-cutting cutter flutes
- ✚ Tips, helix coil, twisted

Preventive Maintenance: Regular functional check and visual inspection. In case of malfunction, the tool must be changed and not used.

#### **Functional Checking For Reamers**

Description and Function: Intramedullary Reamers, Head of the Reamers, Reamer shafts.

Possible Failure Modes:

- ✚ Defective coupling end
- ✚ Dull and non-cutting cutter flutes
- ✚ Tips, helix coil of reamer shaft deformed

Preventive Maintenance: Regular functional checking and visual inspection. In case of malfunction, the tool must be changed and not used.

#### **Functional Check for Screwdriver Blades**

Description and Function: Various designs of driver-link screwdrivers with or without self-tapping function.

Possible Failure Modes:

- ✚ Blade deformation (bent)
- ✚ Deformation of the blade (round)
- ✚ Breaking the non-functional self-guiding mechanism of the blade
- ✚ The corresponding drive connection of the screw head is round or worn.

Preventive Maintenance: Use a suitable tool spray for the mechanism of self-tapping screwdrivers. Regular functional check and visual inspection. Replace in case of failure and do not use. Normal functional check of the corresponding items (screwdrivers - screws).

#### **Functional Check for Torque Limiters**

Description and Function: Includes all torque limiter or display hand tools with or without opening mechanism.

Possible Failure Modes:

- ✚ Malfunction due to wear
- ✚ Corrosion or contamination

Preventive Maintenance: Use a suitable tool spray for the mechanism of self-tapping screwdrivers. After use or cleaning, check all voids for debris and remove them if necessary. In case of malfunction, the tool must be changed and not used.

#### **Functional Check for Scissors and Cutting Tools**

Description and Function: All tools used for cutting tissue and bone material (such as scissors or osteotomes), as well as all tools (such as cutters), are used to cut implants.

Possible Failure Modes:

- ✚ The cutting edge is damaged, the cutting function is no longer fully provided.
- ✚ Spring damaged, articulation not working properly

Preventive Maintenance: Use a suitable instrument spray for the mechanism of all moving parts and articulated surfaces. Handle tools with care. The cutting edge of the tool should be examined before clinical use. The cutting function can be restricted if the cutting edge or tip of the tool is damaged. It recommends the use of silicone-free, mineral-free oil-based lubricant for the maintenance of articulated tools. In case of malfunction, the tool must be

changed and not used.

### **Functional Check for Targeting Handpieces**

Description and Function: Target Handpieces targeting the locking holes of implants (nails, plates)

Possible Failure Modes:

- ✚ Signs caused by a crash of the device
- ✚ Cracks in the polymer
- ✚ Gear damage
- ✚ Deformation of nail adapter
- ✚ Correction pin displacement

Preventive Maintenance: Use a suitable instrument spray for the mechanism of all moving parts and articulated surfaces. Handle with care. Do not hit target hand tools. In case of malfunction, the tool must be changed and not used.

### **Functional Check for Drill Guides**

Description and Function: Soft tissue-protective covers used during drilling

Possible Failure Modes:

- ✚ The scratched outer surface notches on the ends of the sheath

Preventive Maintenance: Use a suitable instrument spray for the mechanism of all moving parts and articulated surfaces. In case of malfunction, the tool must be changed and not used.

### **Functional Check for Forceps, Clamps, and Holding Instruments**

Description and Function: Repositioning bone fragments; Clamping the wire after the wire is put under tension with the wire tensioner

Possible Failure Modes:

- ✚ Deformed functional surfaces (e.g. teeth and locking mechanism)
- ✚ Span between arms
- ✚ Corrosion between neighboring surfaces
- ✚ Corrosion between matching of different components
- ✚ Corrosion on frequently used functional surfaces
- ✚ Corrosion in laser engraving
- ✚ Corrosion of the clamping disc when stressed can cause the disc to break.

Preventive Maintenance: Use a suitable instrument spray for the mechanism of all moving parts and articulated surfaces. In case of malfunction, the tool must be changed and not used.

Always inspect the wire post and especially the clamping disc before reuse. In case of signs of corrosion or cracks, no wire supports are used.