
## TİBİA ÇİVİLERİ - KULLANIM KILAVUZU

### TIBIA NAILS - INSTRUCTIONS FOR USE

## TÜRKÇE

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## ENGLISH

### Symbols and Their Explanations

-  Read the user manual before use!
-  Warning!
-  It is not sterile!
-  Do not reuse !
-  Store in a dry place!
-  Protect from sunlight!
-  Humidity limitation symbol: Store in environments not exceeding 0-80% humidity!
-  Manufacturer
-  CE marking according to (EU) 2017/745 Medical Device Directive

### References

**(EU) 2017/745** Medical Device Regulation

**EN ISO 20417:2021** Medical devices - Information to be provided by the manufacturer

**EN ISO 15223-1** Medical devices - Symbols to be used in the information to be submitted by the manufacturer - Part 1: General requirements

**EN ISO 17664-1** Processing of health care products - Information provided by medical device manufacturers for the processing of medical devices - Part 1: Critical and semi-critical medical devices

**EN ISO 17665-1** Sterilization of healthcare products - Moist heat - Part 1 - Requirements for the development, validation and routine control of a sterilization process for medical devices

**EN ISO 17665-2** Sterilization of health care products - Moist heat – Part 2 - Requirements for the development, validation and routine control of a sterilization process for medical devices

**EN ISO 14630:2012** Non-active surgical implants - General specifications

**ISO 11607-1:2019** (Packaging for Terminally Sterilized Medical Devices - Part 1: Requirements for Materials, Sterile Barrier Systems and Packaging Systems)

### Product Description

Tibia Intramedullary Nail & Accessories should only be used for bone fixation purposes by specialist healthcare professionals (surgeons) with the necessary application qualifications and full experience and certification in this field. These medical devices and accessories, treatment sets and instruments are manufactured for medical interventions for the treatment of patients and cannot be used for other purposes. Otherwise, manufacturers and distributors cannot be held responsible.

Ortimplant Tibia Intramedullary Nails & Accessories are sold non-sterile and all devices are disposable.

Tibia Intramedullary nails and accessories are manufactured from certified ISO 5832-1 stainless steel and certified ISO 5832-3 titanium alloy raw materials. Titanium and stainless steel (SS) are completely biocompatible, corrosion resistant and non-toxic in the biological environment.

Ortimplant implants are produced from materials that provide the above features for special use. Do not use titanium alloy,

stainless steel and nail components in the same place. Due to the risk of metal fatigue and infection, never reuse Tibia Intramedullary nail and accessory implants.

**Note:** The size, type and dimensions of the Tibia Intramedullary nails and accessories to be used, "application sets" suitable for the selected implant and accessories should be determined by the surgeon who will perform the application according to the bone area on which he will use, the condition of the disease, the gender of the patient and the age of the patient.

Before using Ortimplant Tibia Intramedullary nails and accessories, surgeons and the supporting team should carefully read the instructions for use. All necessary documents and materials will be provided to you by the authorized distributor. The surgeon should be aware of the risks that may arise from the use of these products and inform the patient about these risks as much as possible.

#### **Tibia Intramedullary nail implants and accessories;**

- ✓ Tibia Intramedullary nail, end-cap and 4.5mm Tibia Intramedullary nail locking screw,
- ✓ ANT Tibia Intramedullary nail, end-cap & 4.5mm Tibia Intramedullary nail locking screw,

Specialists are responsible for the consequences of using the wrong implant, not meeting the safety criteria in this user guide, or using inappropriate surgical techniques. In such cases, Ortimplant and its authorized distributor are not responsible to the patient and surgeon.

Tibia Intramedullary nails and accessories produced by Ortimplant can only be used with Ortimplant application sets. The surgeon performing the application; The application sets should strictly avoid the use of other trademarks. When performing such an application, the surgeon should definitely avoid applications where he is not sure that the Tibia Intramedullary nail and accessories, hand tools and tools will be sufficient to ensure bone alignment. Otherwise, the surgeon is responsible for the surgery and postoperative complications. Ortimplant and its distributors are not responsible to patients and surgeons under these circumstances.

During the preparation phase before the surgery, the surgeon must ensure that the Tibia Intramedullary nail implants, application sets and the implants he will use are compatible with each other. Ortimplant and its distributors are not responsible for the negative consequences of surgical application in case of any complications.

Tibia Intramedullary Nail implants produced by Ortimplant are produced in different shapes and sizes according to bone structures, bone anatomical shapes and fracture configurations. It is compatible with Ortimplant product structures, and each product structure is designed according to the specific bone problem.

**Note:** Tibia Intramedullary nails are manufactured from ISO 5832-1 stainless steel and ISO 5832-3 titanium alloy certified raw materials. Do not use titanium alloy and stainless steel plates in the same place.

**Note:** To get the best result; Do not use Ortimplant Tibia Intramedullary nails together with parts produced by other companies. Do not reuse bone plate implants due to metal fatigue and risk of infection.

**Note to Surgeon:** Even if the physician is a knowledgeable intermediary between the patient and the company, the important medical information provided in this document must be disclosed to the patient.

## **Intended Use**

Tibia Intramedullary Nail implants and components are designed and manufactured to be used for fracture fixation in the treatment of adult fractures. It can also be used to stabilize corrective osteotomies and reconstruction procedures and to treat pseudoarthrosis. The main purpose of using such implants is to provide compression and solid fracture fixation in the fracture area. As a result of early mobilization of the patient, the chances of pulmonary, gastrointestinal and cardiovascular disorders are reduced.

## Indication Information

Tibia Intramedullary Nail implants and components are designed and applied for rigid bone fixation on bone fractures, pseudoarthroses and osteotomies.

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**Note:** The size, type and dimensions of the implants and accessories to be used and the "application sets" suitable for the selected implant and accessories should be determined by the surgeon who will perform the application according to the bone area to be used, the structure of the bone, the condition of the disease, the gender of the patient and the age of the patient.

Planning the operation should be based on detailed evaluation of the patient's x-rays. In this way, the appropriate Tibia Intramedullary nail size and final combinations determined by the surgeon should be calculated and determined. The patient's allergic sensitivity to the implant material should be determined. Implants should be available in different sizes. Before placing the appropriate implant, the surgeon must ensure that the implant size is compared to the label to use the correct size implant. Extra labels inside the box are for reporting and ordering purposes.

**Note:** Tibia Intramedullary nail implants, components and application sets suitable for the implant should be selected by the surgeon and his team before the operation, and their operation should be checked and verified.

Tibia Intramedullary nail, accessories and appropriate tool sets will be selected and their proper functioning will be checked by the surgeon and the operation team before the operation. Ortimplant Tibia Intramedullary nail and its components are selected by the surgeon and his team before the operation. The surgeon and the team are responsible for selecting the appropriate Tibia Intramedullary nails and accessories for the appropriate patient. Ortimplant and its distributors cannot be held responsible for any disruptions that may arise from improper implant selection and application.

## Contraindication Information

The following situations include, but are not limited to, contraindications;

- ✓ Development of active infection or risk of infection (immunodeficiency),
- ✓ Indicators of local inflammation,
- ✓ Increased fever and white blood cell mass,
- ✓ Morbid obesity,
- ✓ Pregnancy (until the patient's suitability for surgery conditions is confirmed by experts),
- ✓ mental illnesses,
- ✓ Patients with serious congenital anatomical disorders,
- ✓ The presence of another disease, associated use of a drug or device, may potentially reduce or prevent the function of nail implants (e.g., streptococcus or progressive streptococcal infections, congenital functional abnormality, unexplained semimetal elevation, other disabling conditions).
- ✓ Rapidly progressive joint disease, bone erosion, osteopenia and/or osteoporosis. Osteoporosis and osteopenia are contraindications that may limit the degree of healing due to decreased stabilization and/or mechanical fixation.

- ✓ Cases of metal allergy to titanium, stainless steel,
- ✓ Situations where bone graft or fusion is unnecessary,
- ✓ In cases where the implant components selected for use are too large or too small for the bone area to be applied,
- ✓ Patients who do not have enough soft tissue to cover the operation area, who do not have sufficient bone stock, or who do not have sufficient bone quality,
- ✓ Patients expected to have low physiological performance,
- ✓ Patients with any anomalies in the bone anatomy where the nails are aimed at the implant,
- ✓ Some mental or neuromuscular diseases that may cause complications or fixation failure in the postoperative period,
- ✓ Patients who are not willing to follow postoperative instructions
- ✓ Using the products in situations not specified in the indications,

 **Note:** Due to the high risk of unexpected metal corrosion inside the body, two different metal components should not be used together.

## Potential adverse events and complications

All known scientific safety warnings regarding bone surgery will be taken into account by the surgeon and the team during the implantation of the Tibia Intramedullary nail and accessories. In many cases, negative results may be clinic-related rather than device-related. Potential safety warnings, potential adverse events, and complications related to implantation of Tibia Intramedullary nails and accessories are listed below, but are not limited to;

- ✓ Early or late loosening of screw components in Tibia Intramedullary and accessories due to poor bone quality of the patient or incorrect application,
- ✓ Possibility of disassembly, bending and/or breakage of all components within the body.
- ✓ Development of foreign body reaction (allergy) to implants, formation of rash, corrosion of implants (cracks, friction and / or general corrosion),
- ✓ In patients without adequate soft tissue cover, skin penetration due to pressure on the operation side, increase in fibrous tissue, development of skin necrosis and/or pain, and inflammation of surrounding tissues and bone may occur. Improper use of instruments or incorrect placement of the implant may cause damage to surrounding tissue and nerves.
- ✓ After the surgery, changes in bone smoothness, limb weakening and/or length loss in the relevant body part should be monitored by the surgeon and the patient.
- ✓ Development of local and/or general infection,
- ✓ Paralysis (complete or partial), dysesthesia, hypersensitivity, numbness, paresthesia, appearance of radiculopathy and/or pain, depersonalization, loss of neurological function such as neuroma development, loss of consciousness, tingling sensation and/or visible defects,
- ✓ Disorders of bladder fullness, urinary retention and other urinary system disorders,
- ✓ Scar formation can cause nervous system problems or pressure and pain around the nerve,
- ✓ Fracture, microfracture, resorption, development of bone defect or penetration of bone,
- ✓ Pseudo arthrosis, late and poor union of the fractured bone where the surgical operation was performed,
- ✓ Increased/decreased mobility and function in the adjacent joints of the operated bone,
- ✓ The patient cannot continue his daily activities,
- ✓ Loss or decrease in bone density (possibly a result of a decrease in the load on the bone),
- ✓ Complications of the grafted area, including pain, fracture or wound healing problems,
- ✓ Intestinal obstruction or loss of control, gastritis or other gastrointestinal tract disorders
- ✓ Bleeding, hematoma, seroma, edema, hypertension, vascular occlusion, stroke, excessive bleeding, phlebitis, tissue decay, wound dehiscence, vascular injuries or other cardiovascular system disorders,

- ✓ Reproductive system disorders, sexual dysfunction, progression of respiratory problems (e.g. pulmonary embolism, pneumonia, bronchitis),
- ✓ It has been reported that during implantation of nails, especially around joints, inappropriately sized nails can penetrate the joints and cause chondral lesions, which should be avoided to prevent chondral lesions and prevent early joint arthrosis,
- ✓ Joint surfaces may be damaged as a result of placement of inappropriate products, wrong size nails, screws and guide wires during surgery,
- ✓ Incidence of bone nail and screw implants in the joint area during improper application,
- ✓ Patients should stay away from strong magnetic fields.
- ✓ The patient should immediately consult his physician if he hears abnormal sounds, pain or increased temperature in the surgery area.

**Note:** Additional surgery may be required to correct these potential drawbacks.

## Intended User Profile

Ortimplant Tibia Nails should only be used by healthcare professionals.

Surgical procedures should only be performed by medical specialists or medical residents under the responsibility of specialists who have received adequate training in the anatomical structure or formation where the surgery will be performed and are familiar with surgical techniques.

Before performing a surgical procedure, surgical literature on techniques, complications, and hazards should be reviewed. You should carefully read all operating instructions regarding product safety features before use. All necessary materials and documents must be provided by the authorized distributor.

## General Safety Warnings

- ✓ Read the User Manual before use and take into account the warnings.
- ✓ The products are NOT STERILE! Before first use, sterilize according to the instructions specified in this manual.
- ✓ Do not use it for purposes other than its intended use.
- ✓ Tibia Intramedullary nail implant should not be used for purposes other than those stated above.
- ✓ Choosing the appropriate device, appropriate accessory and appropriate application tool set for the patient is very important. Improper use or misalignment of components may cause loss, bending, breakage, or cracking of the device or bones.
- ✓ Until radiographic confirmation of adequate bone development; Other internal or external support devices should be used (orthosis, plaster, etc.) to stabilize the fracture and minimize the load on the fractured area.
- ✓ High pressure occurs in proximal/distal fractures and osteotomies of long bones. This can cause bending or breakage of implant components due to uneven muscle force, which can negatively impact healing.
- ✓ Do not mix incompatible metals.
- ✓ The device is made of metal materials and does not contain waste materials.

## Precautions

The implant should not be reused. Loading too early may cause devices to stop working. In order to prevent the implant from causing discomfort during use, check the implants and application sets (in terms of suitability and full functionality) before use.

## Operation

### Patient and Implant Selection:

Only patients with indications and no absolute contraindications should undergo surgical procedure. However, the final decision and responsibility belongs to the surgeon who will perform the surgery.

Choosing the appropriate size, shape and design for each patient is important for the procedure to be successful. Metallic surgical implants are subjected to repeated loads during use and are limited in shape, size and design. Unless care is taken in selecting the appropriate implant for the patient, unless the implant is placed properly, and if no precautions are taken to reduce the stress on the implant, such stresses can result in buckling and device loosening, metal fatigue, and eventual fracture before healing. These situations may result in device damage or premature removal.

### Before Operation:

- ✓ Only patients who meet the criteria specified in the indications should be selected.
- ✓ The use of Tibia Intramedullary nails should be avoided in patients who have or are prone to the conditions listed in the contraindications.
- ✓ Care should be taken when transporting and storing Tibia Intramedullary nails and accessories. Nails should not be scratched or damaged. They must be protected against aggressive environmental conditions during storage. Implant packages should be checked carefully for damage before use.
- ✓ Devices of appropriate size, type and number should be used during surgery.
- ✓ Surgeons must become familiar with the various components of the Tibia Intramedullary nails before using the surgical equipment and set up the system to verify that all parts and necessary tools are present before the operation begins. Bone nail implants and accessories (as described in the description) cannot be completed with components produced by other manufacturers. Different types of metal should never be used together.
- ✓ Tibia Intramedullary nails and accessories must be clean and sterile before use. Additional sterile components must be available to prevent undesirable situations.

### During the operation:

- ✓ Extreme attention should be paid to the areas through which the spinal column, nerve roots, and peripheral nervous system pass. Damage to nerves leads to loss of neurological functions. During the surgical procedure, maximum attention should be paid to intracranial, intrathoracic, intrapelvic organs, arterial and venous areas, and these structures should be protected. Arterial or venous injuries can result in death.
- ✓ Tibia Breakage, slippage, or misuse of intramedullary nail components or surgical set may result in injury to the patient or surgical personnel.
- ✓ Bone nail implants should not be repeatedly bent excessively. Scratches or nicks on the implant may reduce its functional integrity, so take due care to avoid such situations.
- ✓ Whenever possible or necessary, an imaging system should be used to facilitate the surgical procedure.
- ✓ Do not use bone bonding agents or solvents as the implant components will be difficult to separate from each other and their safety and effectiveness have not been proven. Heat accumulation occurring during surgery can cause neurological injuries and local bone tissue death.
- ✓ Before suturing the soft tissue, change the Tibia Intramedullary nail position if there is an inappropriate nail implant position.

**After the operation:**

- ✓ 1. The patient should be given instructions on how to use the device. If regional weight distribution is recommended or anticipated, the patient should be informed of possible bending of bone nail implants before bone fusion is completed. If the patient is extremely old, frail, or mobile, the risk of bending, loosening, and impact of the temporary fixation device may increase during postoperative rehabilitation. The patient should be warned to avoid falls or sudden jolts in any position.
- ✓ 2. To achieve maximum surgical success, neither the patient nor the device should be exposed to mechanical vibration or shock that could cause the device to loosen. The patient should be informed about this possibility and, in particular, his/her participation in any sports activities should be restricted. The patient should be advised to be accompanied by those nearby and to use auxiliary support devices (crutches, walkers, orthoses, etc.) while standing up. The patient should especially be advised to avoid turning movements. The patient should be warned not to drink or consume a product containing nicotine, alcohol, or steroids during the bone healing process.
- ✓ 3. The patient should be warned about the loss of movement in the surgical area and informed about how to compensate for the permanent physical limitation in body movement.
- ✓ 4. Nonunion or delayed union of the fracture; creates excessive and repetitive stress on the implant that stabilizes the fracture. Excessive load may cause the metal device to bend, loosen, or collapse due to fatigue. It is important not to move the surgical area until appropriate bone union is achieved by x-ray examination. If the device is loose, bent and/or warped, it must be replaced or removed without serious injury. The patient should be warned of these dangers and monitored closely until the bone consolidates.
- ✓ 5. Especially in high-risk patients, prophylactic antibiotic use is precautionary before subsequent surgical procedures (such as dental surgery).
- ✓ 6. Implant fixation devices are designed to provide stability at the fracture site during normal fracture healing. Once the bone is fused, these devices have no functional effect and can be removed. The final decision to remove the device is up to the surgeon and the patient. Many patients may experience one or more of the following complications:
  - Disruption of implant position (leading to possible injuries),
  - Extra risk of injury due to trauma to the surgical area after surgery,
  - Bending, loosening or breaking of the implant in the surgical area, which may make removal difficult,
  - Pain, discomfort or abnormal movement while the implant is in the body,
  - High risk of possible infection,
  - Bone loss caused by excessive pressure on the bone

**Removing the Implant:**

Tibia Intramedullary nails must be removed by surgeons when there is a problem such as pain or infection. This procedure is performed with anesthesia and also requires a surgical procedure. This procedure is described in the "Surgical Technique of Ortimplant Tibia Intramedullary Nails" guide. The extraction adapter provided by the relevant manufacturer is attached to the tip of the Tibia Intramedullary nail and the nail is removed with the recommended extraction technique.

**Note:** Ortimplant recommends removing the implant after the fracture has healed to prevent refracture of the bone and other complications. Large Bone Plates and their components should not be reused. It is not possible to use a device from one patient to another and should not even be considered.

## Surgical Technique

A separate surgical technique guide is provided to explain the Tibia Intramedullary Nail procedure. The surgeon should be familiar with the surgical procedure before using the products. It is also the surgeon's responsibility to consult with colleagues experienced in the procedures and to become familiar with relevant publications before use. Please review the surgical technique document provided by Ortimeplant for the correct application method.

## Sterilization

Steam Sterilization Method has been validated and is recommended as the sterilization method for our products.

- ! The steam sterilizer must comply with the EN 285 standard and be validated according to the EN ISO 17665 standard.
- ! It must be ensured that the maximum allowable loading is not exceeded according to the sterilizer manufacturer's instructions.
- ! The properties of the water to be used in the sterilization process must comply with EN 285 and ISO 17665 standards.

Sterilization Method	Temperature	Time
Steam Sterilization	134 °C	5 min

## Packaging and Storage Conditions

All Ortimeplant Tibia Intramedullary Nails are cleaned and disinfected according to standards-compliant methods. Cleaning and decontamination involves rinsing with neutral cleaning agents followed by washing with distilled water.

Afterwards, all products are packaged in the clean room.

**Note:** It should be noted that the packages containing Ortimeplant Tibia Intramedullary Nails are not damaged. If a delivery note or invoice is used, check that all products are complete, including application sets, before use. Damaged packages and products should not be used and should be returned to Ortimeplant or the distributor.

Tibia Intramedullary Nails Plates should be stored in an environment with humidity below 80% and the environment should not contain corrosive gas.

## Disposal Conditions

Disposal of the products must be carried out in accordance with the national legislation of the country where the products are used.

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**Manufacturer:** Ortimplant Ortopedi ve Tıbbi Aletler Tic. Ltd. Şti.