

# ORTHOPEDIC INSTRUMENT CONTAINERS and trays


" We produce for your health... "



## INSTRUMENT CONTAINER CATALOG

"We produce for your health..."

Ortimplant provides a complete platform of aluminum perforated boxes, screw racks, silicone holders, surgical instrument trays are produced for all surgical instrument sets, which allow the orthopedic surgical instruments to be stored, transported and sterilized by wrapping surgical textile



**Small Fragment Instrument Sets**

**Large Fragment Instrument Sets**

**DHS/DCS Instrument Sets**



# ORTIMPLANT PRODUCT GROUPS

## ● STERILIZATION CONTAINERS

We produce aluminum sterilization containers suitable for all surgical sets which are convenient, hygienic, safe & which can be used without wrapping with surgical textiles and which can maintain sterilization safety for a long time

## ● TRAUMA INSTRUMENT SETS

Orthopedic surgery instrument sets, including carefully designed surgical instruments aluminum boxes, silicone holders suitable for the placement of the surgical instruments for all orthopedic surgeries, provide great convenience to surgeons during surgery

## ● ORTHOPEDIC STERILIZATION CONTAINERS

Aluminum perforated boxes, screw racks, silicone holders, surgical instrument trays are produced for all surgical instrument sets, which allow the orthopedic surgical instruments to be stored, transported and sterilized by wrapping surgical textile

"With its obsessive Quality policy...."

With its obsessive quality policy and with its open structure to development. Ortimeplant goes on its constitution under today's modern institutional structure



## Screw Rack

## Cannulated Screw Instrument Sets



Our expert team has years of experience in development of business and strategy..

From raw materials to services after sales, Ortimeplant considers to apply Total Quality Management System, With respect to this policy

**Ortimeplant has been certified its Quality with the certificates of CE, ISO 9001, ISO 13485.**



Ortimplant manufactures orthopaedics implants, sterilization containers and accessories in line with world norms with its high technologic workshop.

Ortimplant increases its success in orthopaedic field everyday by manufacturing new and developed products as a result of investigation and research studies of the specialist physicians and experienced engineers according to needs of global orthopaedic market.

Together with modern background and educated technical team Ortimplant aims qualified manufacturing as a goal.

From raw materials to services after sales Ortimplant considers to apply Total Quality Management. With respect to this policy Ortimplant has been certified its Quality with the certificates of CE, ISO 9001, ISO 13485.

With the aim of more rapid and continuous satisfaction of customer desires, Ortimplant educates the employees in advanced stages to increase customer satisfaction under all conditions. With periodically given education seminars Ortimplant is in the aim of reaching institutional quality from personal quality and with direction of this goal Ortimplant succeeded the communication between employees by meeting all management and technical teams together in social and cultural activities.

Ortimplant attributes to national and international sectoral fairs to actively present its current and new products to the medical sector in the world.

With its obsessive Quality policy and with its open structure to developments Ortimplant goes on its constution under today's modern instutional structure.



# INDEX



SMALL FRAGMENT LOCKING / NON-LOCKING PLATE INSTRUMENT SET .....	6
LARGE FRAGMENT LOCKING / NON-LOCKING PLATE INSTRUMENT SET .....	7
DHS & DCS INSTRUMENT SET.....	7
2,7/3,5/4,0 SCREW RACK .....	8
4,5/5,0 SCREW RACK .....	8
2,7-3,5-4,0 MM CANNULATED SCREW INSTRUMENT.SET.....	9
4,5-6,5-7,3 MM CANNULATED SCREW INSTRUMENT.SET.....	9
NEW FEMORAL RECONSTRUCTION INTRAMEDULLARY NAIL INSTRUMENT SET .....	10
TIBIAL INTRAMEDULLARY NAIL INSTRUMENT.SET.....	10
3,0/3,5/4,0/5,0/6,5 MM HEADLESS COMPRESSION SCREW INSTRUMENT SET .....	11
CABLE INSTRUMENT SET .....	11
SCREW REMOVAL INSTRUMENT.SET.....	12
MINI FRAGMENT PLATE INSTRUMENT.SET.....	12
2,4 MM MAXILLOFACIAL PLATE INSTRUMENT.SET.....	13
STERNUM RIB PLATE INSTRUMENT.SET.....	14
HUMERAL INTRAMEDULLARY NAIL INSTRUMENT.SET.....	15
FEMUR ANTEGRADE INTRAMEDULLARY NAIL INSTRUMENT SET .....	15
APFNA FEMUR SHORT NAIL INSTRUMENT.SET.....	16
INTRAMEDULLARY NAIL REMOVAL INSTRUMENT.SET.....	16
ANT TIBIAL INTRAMEDULLARY NAIL INSTRUMENT.SET.....	17
SPINAL PEDICLE SCREW SYSTEM INSTRUMENT.SET.....	18
ANTERIOR CERVICAL PLATE SYSTEM INSTRUMENT SET.....	18
AYFIX EXTERNAL FIXATOR INSTRUMENT.SET.....	19
PIN-WIRE INSTRUMENT SET.....	19
FLEXIBLE REAMER SURGERY INSTRUMENT.SET.....	20
HERBERT CANNULATED SCREW INSTRUMENT.SET.....	20
PELVIC RECONSTRUCTION PLATE INSTRUMENT.SET.....	21
ELASTIC NAIL INSTRUMENT.SET.....	21

# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS


SMALL FRAGMENT LOCK. / NON-LOCK. PLATE INSTRUMENT CONTAINER

03200-0130

BOX SIZE / MM  
363X252X148



CODE	PRODUCT
03200-0100	Set
03200-0130	Aluminium Box
03200-0130-01	Box (without lid)
03200-0130-02	Lid
03200-0130-03	Screw Rack
03200-0130-04	Tray-1
03200-0130-05	Tray-2
03200-0130-06	Tray-3





# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS




## LARGE FRAGMENT LOCK./NON-LOCK. PLATE INSTRUMENT SET CONTAINER

03200-0222



BOX SIZE / MM  
435X285X148

CODE	PRODUCT
O3200-0200	Set
O3200-0222	Aluminium Box
O3200-0222-01	BOX (without lid)
O3200-0222-02	Lid
O3200-0222-03	Screw Rack
O3200-0222-04	Tray-1
O3200-0222-05	Tray-2
O3200-0222-06	Tray-3






## DHS & DCS INSTRUMENT CONTAINER

04200-0522



BOX SIZE / MM  
445X259X155

CODE	PRODUCT
O4200-0500	Set
O4200-0522	Aluminium Box
O4200-052-01	Box (without lid)
O4200-052-02	Lid
O4200-0522-03	Screw Rack
O4200-0522-04	Tray-1
O4200-0522-05	Tray-2
O4200-0522-06	Tray-3



# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS

## 2,7/3,5/4,0 SCREW RACK

OSR0002

BOX SIZE / MM  
286X156X123



CODE	Bottom Perforated
OSR0002	Set



## 4,5/5,0 SCREW RACK

OSR0003

BOX SIZE / MM  
286X156X123



CODE	Bottom Perforated
OSR0003	Set





# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS

## 2,7-3,5-4,0 MM CANNULATED SCREW INSTRUMENT CONTAINER

04200-0612

BOX SIZE / MM  
245X250X152



CODE	PRODUCT
O4200-0600	Set
O4200-0612	Aluminium Box
O4200-0612-01	BOX (without lid)
O4200-0612-02	Lid
O4200-0612-03	Screw Rack
O4200-0612-04	Tray-1
O4200-0612-05	Tray-2
O4200-0612-06	Tray-3



## 4,5-6,5-7,3 MM CANNULATED SCREW INSTRUMENT CONTAINER

04200-7329

BOX SIZE / MM  
421X296X153



CODE	PRODUCT
O4200-7300	Set
O4200-7329	Aluminium Box
O4200-7329-01	BOX (without lid)
O4200-7329-02	Lid
O4200-7329-03	Screw Rack
O4200-7329-04	Tray-1
O4200-7329-05	Tray-2
O4200-7329-06	Tray-3



# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS




## NEW FEMORAL RECONS. INTRAMEDULLARY NAIL INSTRUMENT CONTAINER

01200-0847

BOX SIZE / MM  
500X299X165



CODE	PRODUCT
01200-0800	Set
01200-0847	Aluminium Box
01200-0847-01	Box (without lid)
01200-0847-02	Lid
01200-0847-03	Tray-1
01200-0847-04	Tray-2






## TIBIAL INTRAMEDULLARY NAIL INSTRUMENT CONTAINER

01200-0236

BOX SIZE / MM  
498X300X195



CODE	PRODUCT
01200-0200	Set
01200-0236	Aluminium Box
01200-0236	Box (without lid)
01200-0236-02	Lid
01200-0236-03	Tray-1
01200-0236-04	Tray-2







# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS

## 3,0/3,5/4,0/5,0/6,5 MM HEADLESS COMP. SCREW INSTRUMENT CONTAINER 04200-1233

BOX SIZE / MM  
298X213X124




CODE	PRODUCT
O4200-1200	Set
O4200-1233	Aluminium Box
O4200-1233-01	Box (without lid)
O4200-1233-02	Lid
O4200-1233-03	Tray-1
O4200-1233-04	Tray-2





## CABLE INSTRUMENT CONTAINER

04200-1706

BOX SIZE / MM  
427X220X165



CODE	PRODUCT
O4200-1700	Set
O4200-1706	Aluminium Box
O4200-1706-01	Box (without lid)
O4200-1706-02	Lid
O4200-1706-03	Tray-1
O4200-1706-04	Tray-2






# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS

## SCREW REMOVAL INSTRUMENT CONTAINER

04200-1049

BOX SIZE / MM  
**380X226X117**




CODE	PRODUCT
O4200-1000	Set
O4200-1049	Aluminium Box
O4200-1049-01	Box (without lid)
O4200-1049-02	Lid
O4200-1049-03	Tray-1
O4200-1049-04	Tray-2

## MINI FRAGMENT PLATE INSTRUMENT CONTAINER

04200-0427

BOX SIZE / MM  
**402X215X97**





CODE	PRODUCT
O4200-0400	Set
O4200-0427	Aluminium Box
O4200-0427-01	Box (without lid)
O4200-0427-02	Lid
O4200-0427-03	Tray-1
O4200-0427-07	Mini Box-1
O4200-0427-08	Mini Box-2
O4200-0427-09	Mini Box-3

# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS


## 2,4 MM MAXILLOFACIAL PLATE INSTRUMENT CONTAINER

O3128-0215

BOX SIZE / MM  
**400X250X62**



CODE	PRODUCT
O3128-0200	Set
O3128-0215	Aluminium Box
O3128-0215-01	BOX (without lid)
O3128-0215-02	Lid
O3128-0215-03	Tray-1
O3128-0215-04	Mini Box



BOX SIZE / MM  
353X243X90



CODE	PRODUCT
O4200-0900	Set
O4200-0909	Aluminium Box
O4200-0909-01	Box (without lid)
O4200-0909-02	Lid
O4200-0909-03	Tray-1
O4200-0909-04	Tray-2





# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS


## HUMERAL INTRAMEDULLARY NAIL INSTRUMENT CONTAINER

01200-0131

BOX SIZE / MM  
442X270X165



CODE	PRODUCT
01200-0100	Set
01200-0131	Aluminium Box
01200-0131-02	BOX (without lid)
01200-0131-02	Lid
01200-0131-03	Tray-1




## FEMUR ANTEGRADE INTRAMEDULLARY NAIL INSTRUMENT CONTAINER

01200-0334

BOX SIZE / MM  
510X276X165



CODE	PRODUCT
01200-0300	Set
01200-0334	Aluminium Box
01200-0334-01	BOX (without lid)
01200-0334-02	Lid
01200-0334-03	Tray-1
01200-0334-04	Tray-2
01200-0334-05	Tray-3




# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS

## APFNA FEMUR SHORT NAIL INSTRUMENT CONTAINER

01200-0762

BOX SIZE / MM  
556X306X205



CODE	PRODUCT
O1200-0700	Set
O1200-0762	Aluminium Box
O1200-0762-01	Box (without lid)
O1200-0762-02	Lid
O1200-0762-03	Tray-1
O1200-0762-04	Tray-2

## INTRAMEDULLARY NAIL REMOVAL INSTRUMENT CONTAINER

01200-1225

BOX SIZE / MM  
435X224X115



CODE	PRODUCT
O1200-1200	Set
O1200-1225	Aluminium Box
O1200-1225-01	Box (without lid)
O1200-1225-02	Lid
O1200-1225-03	Tray-1
O1200-1225-04	Tray-2
O1200-1225-05	Tray-3



# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS


## ANT TIBIAL INTRAMEDULLARY NAIL INSTRUMENT CONTAINER

01200-0951

BOX SIZE / MM  
**500X304X164**



CODE	PRODUCT
01200-0900	Set
01200-0951	Aluminium Box
01200-0951-01	Box (without lid)
01200-0951-02	Lid
01200-0951-03	Tray-1
01200-0951-04	Tray-2
01200-0951-05	Tray-3






# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS

## SPINAL PEDICLE SCREW SYSTEM INSTRUMENT CONTAINER

O2200-0131

BOX SIZE / MM  
**552X274X175**




CODE	PRODUCT
O2200-0100	Set
O2200-0131	Aluminium Box
O2200-0131-01	Box (without lid)
O2200-0131-02	Lid
O2200-0131-02	Tray-1
O2200-0131-03	Tray-2
O2200-0131-04	Tray-3

## ANTERIOR CERVICAL PLATE SYSTEM INSTRUMENT CONTAINER

O2200-0223

BOX SIZE / MM  
**353X243X150**



CODE	PRODUCT
O2200-0200	Set
O2200-0223	Aluminium Box
O2200-0223-01	Box (without lid)
O2200-0223-02	Lid
O2200-0223-02	Tray-1
O2200-0223-03	Tray-2

# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS


## AYFIX EXTERNAL FIXATOR INSTRUMENT CONTAINER

05200-0122

BOX SIZE / MM  
411X259X120



CODE	PRODUCT
O5200-0100	Set
O5200-0122	Aluminium Box
O5200-0122-01	Box (without lid)
O5200-0122-02	Lid
O5200-0122-03	Tray-1
O5200-0122-04	Tray-2



## PIN-WIRE INSTRUMENT CONTAINER

06200-0107

BOX SIZE / MM  
396X251X115



CODE	PRODUCT
O6200-0100	Set
O6200-0107	Aluminium Box
O6200-0107-01	Box (without lid)
O6200-0107-02	Lid
O6200-0107-03	Tray-1
O6200-0107-04	Tray-2







# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS

## FLEXIBLE REAMER SURGERY INSTRUMENT CONTAINER

O5201-0116

BOX SIZE / MM  
500X165X62





CODE	PRODUCT
O5201-0100	Set
O5201-0116	Aluminium Box
O5201-0116-01	BOX (without lid)
O5201-0116-02	Lid

## HERBERT CANNULATED SCREW INSTRUMENT CONTAINER

O4200-1509

BOX SIZE / MM  
292X217X75



CODE	PRODUCT
O4200-1500	Set
O4200-1509	Aluminium Box
O4200-1509-01	Box (without lid)
O4200-1509-02	Lid
O4200-1509-03	Tray-1



# ORTHOPEDIC INSTRUMENT CONTAINERS & TRAYS


## PELVIC RECONSTRUCTION PLATE INSTRUMENT CONTAINER

04200-0831

BOX SIZE / MM  
537X315X165



CODE	PRODUCT
O4200-0800	Set
O4200-0831	Aluminium Box
O4200-0831-02	Box (without lid)
O4200-0831-02	Lid
O4200-0831-03	Tray-1
O4200-0831-04	Tray-2




## ELASTIC NAIL INSTRUMENT CONTAINER

01200-0618


BOX SIZE / MM  
504X352X125




CODE	PRODUCT
O1200-0600	Set
O1200-0618	Aluminium Box
O1200-0618-01	Box (without lid)
O1200-0618-02	Lid
O1200-0618-03	Tray-1



# PRODUCT GROUPS



TRAUMA PRODUCTS



TRAUMA INSTRUMENT SETS



ORTHOPEDIC SURGICAL INSTRUMENTS



STERILIZATION CONTAINERS



STERILIZATION BASKETS




ORTHOPEDIC STERILIZATION CONTAINERS





**Ortimplant**  
**manufactures orthopaedics implants,**  
sterilization containers and accessories in line with world norms  
with its high technologic workshop.





# ortimplant

---

[www.ortimplant.com](http://www.ortimplant.com)

Ortimplant Samsun  
Kirazlık Mahallesi, 1034. Cadde No:14 Tekkek y  
SAMSUN / TURKEY  
[export@ortimplant.com](mailto:export@ortimplant.com)

Workshop  
Organize Sanayi Sitesi Vali Mehmet Erdo an Cebeci  
Bulvarı No:34 Tekkek y / SAMSUN / TURKEY

Tel : +90 362 266 97 28  
Fax : +90 362 266 58 72  
[export@ortimplant.com](mailto:export@ortimplant.com)

“  
**Quality  
in  
Health**